

Product datasheet for SC400292

MIR223 Human MicroRNA Expression Plasmid (MI0000300)

Product data:

Product	MicroRNA Expression Plasmids
Type:	
Species:	Human
Vector:	pCMV-MIR (PCMVMIR)
Precursor:	CCUGGCCUCCUGCAGUGCCACGCUCCGUGUAUUUGACAAGCUGAGUUGGACACUCCAUGUGGUAGAGUGUCAGUUUGUCAAAUACCCCAAGUGCGGCACAUGCUU
Mature Sequence:	CGUGUAUUUGACAAGCUGAGUU
Full insert:	>MI0000300 TAGAGGAACCCACTGAAGCTTGGGACCCCTAGAAATGCCATAGCTACAGGACTCAGGCCTGATGAGCTTC CAGCTGAGCACTGGGTGTCACCTCGGGCTTTACCTGCTTATCTTCAGGATCTCTCTTCTGGTTAGGAACCC TATGCCTATTTTGTCACTTCCCTATTCTGGTGCTTTGGTTGGTCCTTTGAGCAAGATCCCCGGGGCTAA GGGTGTGACTTCATCATTCTTTCTCTCTTTCCCTCTAGGGTCACATCTCCAGGAAGATCTCACTTC CCCACAGAAGCTCTTGGCCTGGCCTCCTGCAGTGCCACGCTCCGTGTATTTGACAAGCTGAGTTGGACAC TCCATGTGGTAGAGTGTGAGTTGTCAAATACCCCAAGTGCGGCACATGCTTACCAGCTCTAGGCCAGGG CAGATGGGATATGACGAATGGACTGCCAGCTGGATACAAGGATGCTCACCAAGCACCAAGTTCTCACAAG TTATTTTATGTGACTTTGCAGGAAGTGGGACTTATATCTGAGGACACCAGGGGAAAAGTGTGGCATCTC AGGGAAATACAGCCCTGGGCTGTGTCTACACACCATGAGAGTGCTGATGGGGCGCAATAGTCTTGAA AATGTATAAAGTGTCCAGGAATGGAAGTGCTTTGATTATTATTTTCTCTTCATATTTCCCTC CCAG
OTI	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence imiRBase . The
Disclaimer:	flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0000300
Synonyms:	MIRN223 mir-223 miRNA223



[View online »](#)

Locus ID: 407008