

Product datasheet for **SC400297**

MIR24-2 Human MicroRNA Expression Plasmid (MI0000081)

Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Human
Vector:	pCMV-MIR (PCMV MIR)
Precursor:	CUCUGCCUCCCGUGCCUACUGAGCUGAAACACAGUUGGUUUGUGUACACUGGCUCAGUUCAGCAGGAACAGGG
Mature Sequence:	UGCCUACUGAGCUGAAACACAG
Full insert:	>MI0000081 GGATGGGATTTGCTTCCTGTGCACAAATCACATTGCCAGGGATTCCAACCGACCCTGAGCTCTGCCACCG AGGATGCTGCCCGGGGACGGGGTGGCAGAGAGGCCCGAAGCCTGTGCCTGGCCTGAGGAGCAGGGCTTA GCTGCTTGAGCAGGGTCCACACCAAGTCGTGTTACAGTGGCTAAGTTCCGCCCCAGGCCCTCACC TCCTCTGGCCTTGCCGCCTGTCCCTGTGCGCCTGTCTGCCTGCCATCCTGTGCCTGGCCTCCCTGG GCTCTGCCTCCCGTGCCTACTGAGCTGAAACACAGTTGGTTTTGTGTACTGGCTCAGTTCAGCAGGAAC AGGGGTCAAGCCCCCTTGAGCCTGCAGCCCCTGCCTCCCTGGGTGGGCTGATGCTGGAGCAGAGATG AGGACTCAGAATCAGACCTGTGTCTGGAGGAGGGATGTGGTGGGTGGGTTGGCTGGGCCAAATGTGTG CTGCAGGCCCTGATCCCCAACTGCAACTGGGACCCCTGCATGGCCACAGCTCAGGCTGGGCTGTGGT GCCAGCATAGATAGGTGGGTGAGTGGGTGGCCCTTCCATTAAGGGAAGCCAGCTGTGTCTTTCCG
OTI Disclaimer:	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in miRBase . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	MI0000085 Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0000081
Synonyms:	MIRN24-2 miR-24-2 miRNA24-2
Locus ID:	407013



[View online »](#)