

## Product datasheet for **SC400303**

### MIR27B Human MicroRNA Expression Plasmid (MI0000440)

#### Product data:

**Product Type:** MicroRNA Expression Plasmids

**Species:** Human

**Vector:** pCMV-MIR (PCMVIR)

**Precursor:** ACCUCUCUACAAGGUGCAGAGCUUAGCUGAUUGGUGAACAGUGAUUGGUUCCGCUUUGUUCACAGUGGCUAAGUUCUGCACCUGAAGAGAAGGUG

**Mature Sequence:** AGAGCUUAGCUGAUUGGUGAAC

**Full insert:** >MI0000440 GGAAGCCCAGTGTGTGCAGACAGCACGGGGTGGCGCTGCTCTCAGGTGCTCTGGCTGCTTGGGTTCTCTGG  
CATGCTGATTTGTGACTTAAGATTAATAACATTGCCAGGGATTACCACGCAACCACGACCTTGGCTGC  
TCCTCCAGAAACCGTGGTCGCGTCACTGCAGATTGGAGAACAGGTGCATCTCGTAGCTCTTCTTTGGAA  
ACAAAAGAAGCCACCAGCTGAGGAAGATGCTCACCGGTCACCGTCCCTTTATTTATGCCAGCGATGACC  
TCTCTAACAAAGGTGCAGAGCTTAGCTGATTGGTGAACAGTGATTGGTTTCCGCTTTGTTACAGTGGCTA  
AGTTCTGCACCTGAAGAGAAGGTGAGATGGGGACAGTTAAGTTGGAGCCGCTGGGGCAGAGGCCGTTGCT  
GACGGGCCGCGCTGCTGCACAGTCACTTGGGTGCGGAGCGGATCCTGGAGGATGAGAGACCACTTG  
ACCCCAAGGATGCACTGTCTCCTGCTGGGAATGCTAGCCATGACTGAGTCCCTTAAGTCTGTCCACAGAA  
ACATGCACTAATCGGACATCTGTCTGAAAGGTCAAATGTATTGAAAGTTGCAAAAATTCTTCTTACAAAA  
AACTAAAACCAATGCATCACCTAAGTCGTGTGAAATC

**OTI Disclaimer:** All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

**E. coli Selection:** Kanamycin (25 ug/mL)

**Mammalian Cell Selection:** Neomycin

**Components:** 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

**Also Contains:** Both mature and minor (star) miRNA can be over expressed from the construct.

**RefSeq:** [MI0000440](#)

**Synonyms:** MIRN27B miRNA27B

**Locus ID:** 407019



[View online »](#)