

Product datasheet for **SC400316**

MIR302a Human MicroRNA Expression Plasmid (MI0000738)

Product data:

Product Type: MicroRNA Expression Plasmids

Species: Human

Vector: pCMV-MIR (PCMVIR)

Precursor: CCACCACUUAACGUGGAUGUACUUGCUUUGAAACUAAAGAAGUAAGUGCUUCAUGUUUUGGUGAUGG

Mature Sequence: ACUUAACGUGGAUGUACUUGCU

Full insert: >MI0000738 TTAACATGGAAGTGCTTCTGTGACTTTAAAAGTAAGTGCTCCATGTTTTAGTAGGAGTGAATCCAATT
TACTTCTCCAAAATAGAACACGCTAACCTCATTGAAGGGATCCCTTTGCTTTAACATGGGGGTACCTG
CTGTGTGAAACAAAAGTAAGTGCTTCCATGTTTCAGTGGAGGTGTCTCCAAGCCAGCACACCTTTTGTTA
CAAAATTTTTTTGTTATTGTGTTTTAAGGTTACTAAGCTTGTACAGGTTAAAGGATTCTAACTTTTTCC
AAGACTGGGCTCCCCACCACTTAAACGTGGATGTACTTGCTTTGAACTAAAGAAGTAAGTGCTCCATG
TTTTGGTGATGGTAAGTCTTCTTTTACATTTTATTATTTTTTAGAAAATAACTTTATTGTATTGACC
GCAGCTCATATATTTAAGCTTTATTTTGTATTTTACATCTGTTAAGGGGCCCTCTACTTTAACATGG
AGGCACTTGCTGTGACATGACAAAATAAGTGCTTCCATGTTTGAGTGTGGTGGTTCCTACCTAATCAGC
AATTGCGTTAACGCCACACTGTGTGCAGTTCCTGGCTACAGGCCATTACTGTTGCTAATATGCAACTCT
GTTGAATATAAATTGGAATTGCA

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0000738](#)

Synonyms: MIRN302 MIRN302A hsa-mir-302

Locus ID: 407028



[View online »](#)