

## Product datasheet for SC400332

### MIR320B-2 Human MicroRNA Expression Plasmid (MI0003839)

#### Product data:

**Product** MicroRNA Expression Plasmids

**Type:**

**Species:** Human

**Vector:** pCMV-MIR (PCMVMIR)

**Precursor:** UGUUAUUUUUUGUCUUCUACCUAAGAAUUCUGUCUCUUAGGCUUUCUCUCCAGAUUUCCCAAAGUUGGAAAAGCUGGGUUGAGAGGGCAAAGGAAAAA

**Mature** AAAAGCUGGGUUGAGAGGGCAA

**Sequence:**

**Full insert:** >MI0003839 GTTAAAAATTGAGGGTCCATTATAACAAAATAATGTTTTCTAAAGTTAAAAGTTGCTACTGTTACCTC  
 TCAATTTTAACTTTTTTTTTCTTTTTTAATTAATGTTTTTACCCATGGCAAGCTGTAATAGCTTTTTTG  
 AGGGGAGGTAGGTGCTTTATAAAGAACAGTAGGTGCTGCTTATCAACAGATGAAAGGAGGGTTCTTTTC  
 AGGCAACCATCTCATTTGTGAGTGAATGCACCTTCTCTTTAAAGTGCTGGGATTGTTAGTGCCATTTTAA  
 TTGTAATATCAGAATTGTTATTTTTGTCTTCTACCTAAGAATTCTGTCTCTTAGGCTTTCTCTCCCA  
 GATTTCCCAAAGTTGGGAAAAGCTGGGTTGAGAGGGCAAAGGAAAAAAGAATTCTGTCTCTGACAT  
 AATTAGATAGGAACCAAGTTGGGAAGCTGTAAGAATAATGCAGGTGCAAGGTGGTGGTGGTTGAGCCGG  
 GTGATAGCTGTGGATGTAGAAGGAATCTGAATATATTGTGTCATAGGGTTGACCTGATTTGCTAATGGAG  
 TAGTTAAGGATGTGGGAAAGTGGAATCAAGCATGGCTCCAAGGTCTGGGCCTGAAAACCGGGAGAATGA  
 GTTCACATTAAGTAAAGACGGGAAAGACAATGGTAGGGCCTGTTTAGGGAAGAGCTAGGAAGAGATTAGC  
 ATCTCATTTCTAATGATGTAATACATAG

**OTI** All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence [imiRBase](#). Th  
**Disclaimer:** clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the fun

**E. coli** Kanamycin (25 ug/mL)

**Selection:**

**Mammalian** Neomycin

**Cell**

**Selection:**

**Components:** 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

**Also** Both mature and minor (star) miRNA can be over expressed from the construct.

**Contains:**

**RefSeq:** [MI0003839](#)



[View online »](#)