

Product datasheet for **SC400338**

MIR324 Human MicroRNA Expression Plasmid (MI0000813)

Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Human
Vector:	pCMV-MIR (PCMVMIR)
Precursor:	CUGACUAUGCCUCCCCGAUCCCUAGGGCAUUGGUGUAAAGCUGGAGACCCACUGCCCCAGGUGCUGCUGGGGUUGUAGUC
Mature Sequence:	CGCAUCCCUAGGGCAUUGGUGU
Full insert:	>MI0000813 TAGCACTCACCATGTAAGCCATGGACTGAGGTTGCATAGTTGGGACATGGGAAGGAAAATTGCAAAGGGC TTTGTGCACTTGGCCTCATCACCCAGATCTCCAAGATAAGGGCTGACCTAGCTTGTGCAAGTCAAGCAGCAG TACTTGTCTGGGTCACTTTCATCAGGTGCTTCCAGGATTTGTTTTCTAAAAGGGGTGGATGTAAGGGA TGAGGTAGAATTAACCTCTGGTACTGCTGGCAGGCACCTGAGCAGAACATCATTGCTGTCTCTCTTCGCA GAAGCTGAGCTGACTATGCCTCCCCGCATCCCCTAGGGCATTGGTGTAAAGCTGGAGACCCACTGCCCA GGTGTCTGCTGGGGTTGTAGTCTGACCCGACTGGGAAGAAAGCCCCAGGGCTCCAGGGAGAGGGGCTTG GAGGCCCTCACCTCAGTTACATACTGCAGCATAACCATCCGTGCCAGCTTCTCCTGGATCAGCCCAAAGT TGTGAATTTTCTCCCCAACTGGGTACGATTAGTGGCATGATCTACCTGGAAGAGGGTCCACACATCCCC CTGTGGTTCAGTGTGGTCTGCAGTCTCCCTAGGAGAGGGCTGGGCTTGCGCCAGAGGGATGGGTTTTG CATACAACCAGAGTTCAGTCTACGCTGGATGGAGACCTAAG
OTI Disclaimer:	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in miRBase . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0000813
Synonyms:	MIRN324 hsa-mir-324
Locus ID:	442898



[View online »](#)