

Product datasheet for **SC400360**

MIR362 Human MicroRNA Expression Plasmid (MI0000762)

Product data:

Product MicroRNA Expression Plasmids

Type:

Species: Human

Vector: pCMV-MIR (PCMVMIR)

Precursor: CUUGAAUCCUUGGAACCUAGGUGUGAGUGCUAUUUCAGUGCAACACACCUAUUCAAGGAUUCAAA

Mature Sequence: AAUCCUUGGAACCUAGGUGUGAGU

Full insert: >MI0000762

ATTCTCATCATCCATTACACCTACATGCCCTGAGAGTGGGTGGCCACAGAGTGGACTTCTCATTACCACTGTCACCCACAAACATACACAAAAGCATGAATCT

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence [imiRBase](#). The

E. coli Selection:

Kanamycin (25 ug/mL)

Mammalian Cell Selection:

Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0000762](#)

Synonyms: MIRN362 hsa-mir-362

Locus ID: 574030



[View online »](#)