

Product datasheet for **SC400437**

MIR509-2 Human MicroRNA Expression Plasmid (MI0005530)

Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Human
Vector:	pCMV-MIR (PCMVMIR)
Precursor:	CAUGCUGUGUGGUACCCUACUGCAGACAGUGGCAUCAUGUAUAAUUAUAAAAUGAUUGGUACGUCUGUGGGUAGAGUACUGCAUGACAC
Mature Sequence:	UACUGCAGACAGUGGCAAUCA
Full insert:	>MI0005530 TTATAAGGAACAGGATAAGACATGCAAGGATATTTCTTCTGTGAGTAACAGAACCTAACAGAAATGTTTA ATATTTCTCAGTTCTGAAAACGTGCATGTGAAAGATGCAACGTGGGAGTGGACAGGACTCAAAGCGGGAA ATATAAGGAAAGATAGTGGATTGAATTCTGCAGAAGTCCTTCGTCCAAATTCATGGCCACGCACAGCT TGTAATATTGAATGGTTGTTTTAAAGGTGAATGGTGGGTATTAAGGCAAGGCTGCCATCCTCAGACAT GCTGTGTGTGGTACCCTACTGCAGACAGTGGCAATCATGTATAATTAATAATGATTGGTACGTCTGTGGG TAGAGTACTGCATGACACGTGCAACATACATGATGACTGTGTGTGTGGAGGCATTTAGTTGCATG CAGAGGGATGCTTTCCTTACAGTGCAGGTATTTACCCTCACAACCCTCCACTTCCATGCTGGACTACC TCAGAAAACCTCTCAGGTGCCTGACCACCAGACACGTGATCTAGGAATTTAAGGTCTTGCAACAAGCC TCCATTGAGTTCTTCTGAGATTGATTTCCATTCTATGATGTTGCCTGATTGAACCAAGGTTTATCTTG CCATTTTACAGTTTCTTCTGCCT
OTI Disclaimer:	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in miRBase . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	<u>MI0003196</u> Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	<u>MI0005530</u>


[View online »](#)