

Product datasheet for **SC400442**

MIR512-1 Human MicroRNA Expression Plasmid (MI0003140)

Product data:

Product Type: MicroRNA Expression Plasmids

Species: Human

Vector: pCMV-MIR (PCMVMIR)

Precursor: UCUCAGUCUGUGGCACUCAGCCUUGAGGGCACUUUCUGGUGCCAGAAUGAAAGUGCUGUCAUAGCUGAGGUCCAUGACUGAGG

Mature Sequence: CACUCAGCCUUGAGGGCACUUUC

Full insert: >MI0003140 CTGGGCAACAGAGCAAGACTAAGAAAAAAAAAGAAAGAAAAAGAAAAAGAAAAACAAAGTAATGTAAAAGG
TGTTATTCAGGCTGGGGACCCCAAAGTGCTGGGATTACAGGCATGAGCCCCAGGGCAGGCCTAGGATTT
GCATTACAAGGCTCTGCATTGTTCTGAGAACTGGAGAGAGTGCTCAATCTACCTTTTCGCTATTGAGCAA
CATTACAGAAGTCAGTTTTAAATCTCTTATCTTCTGCATAGAGGATGTGCCTGCAGTTTCCAGGTACTGGA
CAATACAGGGAGGGTACTTCTCAGTCTGTGGCACTCAGCCTTGAGGGCACTTTCTGGTGCCAGAATGAAA
GTGCTGTCATAGCTGAGGTCCAATGACTGAGGCGAGCACCGAAGAAACACCATGGGGGGGAGGGGGGGCG
GGGGGCTCCAGGAGCCACTGCAGGTAAAAGCTAAACCGTGGTGCTTGTCTGCGCCATAAACTGGACTTC
AGCCCAGCTTTAAGGAAATGACCACGGTCTCTGGGTTGTGAGGCTGGGTTCCGTAGCAGCTATAATGCTT
GGATTACAGGGATAGCCTGGGTACACCATTCAGAGAGAAGAGGCAACCAATCCAGACCCCAAAAGCGGGT
TCTTGGATCTCTTGCTGGAAAGAATTCGAGGCAAGTCACAACACAGCGAA

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: [MI0003141](#)
Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0003140](#)



[View online »](#)