

Product datasheet for **SC400445**

MIR513A2 Human MicroRNA Expression Plasmid (MI0003192)

Product data:

Product MicroRNA Expression Plasmids

Type:

Species: Human

Vector: pCMV-MIR (PCMVMIR)

Precursor: GGAUGCCACAUAUCAGCAUUCAGUGGAGGAGGUGUCAUUUAUGUGAACUAAAAUUAUUUACCUUUCUGAGAAGGGUAAUGUACAGC

Mature UUCACAGGGAGGUGUCAU

Sequence:

Full insert: >MI0003192 TGAGGCAGGAGAAGTGCTTAGACACAGAAGGGGGAGGTTGCAGTGAGCCGAGATCACACCAATGCATGCC
AGTCTGGGCGACAGAGCCAGACTTCATCTGAAAACAAAAACAAAAACAAAAATCTGCATGTAAA
ATCTCACTTGGGGCAGGATGAAAAGGGGAGCGCCGGGGGCTACAAGTTGCATTGTCCCTTGGCATAAATC
TTGGTGGCCATGCCAAGGCATACACCTATAGGACTTGAAAAAACGGCATGAATAGGGAGCATTGTGTT
CTTGATGCCACATTCAGCCATTCAAGTGTGCAGTGCCTTTCACAGGGAGGTGTCATTTATGTGAATAAA
ATATAAATTTACCTTTCTGAGAAGGGTAATGTACAGCATGCACTGCATATGTGGTGTCCCTGAGGATGG
GGGTGTAGCTGCGAGCATGGTGTGCTGTTTTGCATCCTGTGCACCCACTCTTCAGTGCTGGGCTACTT
CAGGGAACCTTTTGTGCTTCTAAATGACTATACATAGGAAGTGTGGACTCCTGAGAAAACGCATGCACT
TCGGTTCCTTAAGAATTTATCTATTTTGCCTAATGAGGCTTGGCTCAAATAAGAAATACATTGGATAGA
GGCCACAGCCATTCTCGTCAGCTGGAAGGCTATATTCCAGATTCTAACTTTTCTCACTTA

OTI All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence [imiRBase](#). The
Disclaimer: a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect mature miRNA.

E. coli Kanamycin (25 ug/mL)

Selection:

Mammalian Neomycin

Cell

Selection:

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Both mature and minor (star) miRNA can be over expressed from the construct.

Contains:

RefSeq: [MI0003192](#)



[View online »](#)