

Product datasheet for **SC400446**

MIR513B Human MicroRNA Expression Plasmid (MI0006648)

Product data:

Product Type: MicroRNA Expression Plasmids
Species: Human
Vector: pCMV-MIR (PCMVMIR)
Precursor: GUGUACAGUGCCUUUACAAGGAGGUGUCAUUUAUGUGAACUAAAAUAAAAUGUACCCUUUUUGAGAGGAGUAAUGUACAGCA
Mature Sequence: UUCACAAGGAGGUGUCAUUUAU

Full insert: >MI0006648 TAACCCAGGCTTGAGTATAGTGGCCTGATCTCCTCTCACTGAAACCTGACTCCCAGGGTCAAGGTGTACT
 CCTGCCTCAGCTCCTGGGTGGCTGGGATTATAACAAAATAAAACAAAATCTGCATGTTAAAATCTCACT
 CAGGGCAGGATGAAAGGGGAGCACTGGGGCTAGGATTTGCATAATTCCTTGGCATAAATCTTGGTGGCC
 ATGCTAATGTATACACCTATATAGGGACTTGAATAAAAATGGCATGAATAGGGAGCATTTGGTCTGGGAT
 GCCACATTCAGCCATTCAGTGTACAGTGCCTTTCACAAGGAGGTGTCATTTATGTGAACATAAATATAAA
 TGTACCTTTTGGAGAGGAGTAATGTACAGCATGCACTGCATATGTGGTGTCCCTAGGGATGGGGCTGTA
 GCTGCGAGCATGGTGTGCCTGTTTTGGATCCTGTGCACCCACTCCTCAGTGTGGGCTACTTCAAGGAA
 CTTTTGTGCCTTCTAGACGACTAGACATAGGAAGTGTGGACTCCTGAGAAAACGCATGCACTCAAGTTC
 TTTAAGAATTTAACTATTTTGCCTAAATGAGGCTTGGCTCAAATAAGAAATACATTGTATAGAGGCCACA
 GCCATTCTTCGTCAGCTGGAAGGTTGTACTCC

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: [MI0003191](#)
 Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0006648](#)


[View online »](#)