

Product datasheet for **SC400452**

MIR515-2 Human MicroRNA Expression Plasmid (MI0003147)

Product data:

Product Type: MicroRNA Expression Plasmids
Species: Human
Vector: pCMV-MIR (PCMVMIR)
Precursor: UCUCAGCAGUCAUUCUCCAAAGAAAGCACUUUCUGUUGUCUGAAAGCAGAGUGCCUUCUUUUGGAGCGUACUGUUUGAGA
Mature Sequence: UUCUCCAAAGAAAGCACUUUCUG

Full insert: >MI0003147 GAAATCCAGCCCAATTCTCATGGTTAAATATAAGGTATGTGTGGTAGGCATTGCTTTTTCTCTTTGGGGA
 CAGAACTCAGGAGGATTGTCTTTGATGAACAAGGCTAACCTGCTGAGCCTTTGAAGCAAGGAACTGGAG
 ATGGTCCTTTTAGGGGTTTATATTCTGGATTCCAGAAAACATGCAGAGAGAGTCAATAAATGTATCTTTA
 TTTTGTGTCCATTTTAACCAGGTGAGGAAAATTCCGTCAAAAACCCAGACTTTTGGAGCGAGAAGATC
 TCATGCAGTCATTCTCCAAAAGAAAGCACTTTCTGTTGTCTGAAAGCAGAGTGCCTTCTTTTGGAGCGTT
 ACTGTTTGAGAAAAACCACGTTGAAGTTGATGCTGATCTTGGTAACACATTTGCAGAGCGTGCTTATCAT
 CAGACTTGCATGATGTTGGGGTTCTGTTTTGTTAGTTGTTTGCATCATGGGGTCTCTGTTGCCCAGG
 CTGGAGTGGGTGGCGCTTCCAACCTAGATCTCTTGGGCTCAAGTGGTCCCGTGCTTGTCTTTTTTTTTT
 TTCTTGAAACAGAGTCTCTGTCTCTGCCCAGGCTTGAGCGCAGTGCCATGTTATCTGCTCACTGCAACC
 TCTCCTCCCAAGTAGATGG

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: [MI0003144](#)
 Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0003147](#)


[View online »](#)