

Product datasheet for SC400515

MIR548h-4 Human MicroRNA Expression Plasmid (MI0006414)

Product data:

Product MicroRNA Expression Plasmids

Type:

Species: Human

Vector: pCMV-MIR (PCMVMIR)

Precursor: GCUAUUAGGUUGUGCAAAGUAAUCGCGUUUUUGUCAUUACUUUAAUUACUUUACGUUUAUUAAUGACAAAACCGCAAUUACUUUUGACCAACCUAAUAC

Mature AAAAGUAAUCGCGUUUUUGUC

Sequence:

Full insert: >MI0006414 GCTTGTTCATGAAAGAATAAGACAGTTTAAATTTGAGCTCATTCTGCCTTTTCTTTTAAATAGCT
 TTATCAAAGTATAAATTGATAGACAGAAAAACCCACATGTATTTATGTGTACAATTTGATGAGCTTGG
 GTATACATATACCCCATGATATCATCACAATCAAAATGATAACATATCCATCGCTCCAAAAGTTTCC
 TCTGCCCCCTGCTTTTGTGGCAAGTGACACAGAAAACCTGAAGCCAACTTGATCCTTCGTTGTT
 CCATTTCCAAAGCTATTAGTTGGTGCAAAAGTAATCGCGTTTGTGATTACTTTAATTACTTTACGT
 TTCATTAATGACAAAACCGCAATTACTTTGCACCAACCTAATACTTGCTATCATCATTCTTGTCTAT
 TTTCCATTCTCTTATAAATGTTGTACAACTAAATCACACTTAAATAAATAAATGGGAAATAAGAAA
 AAGGAAAAAAGGATGAATGCAGGAAAAAACAAGTGCCTGGGTCTGTCTTCTGTGATTCCCAAAATAACA
 AGCACAAGCCCTGATATTCATATCCAGCATTACCAAGTTGCCCTTTTTTTTTTGTGATGGAGTTTC
 ACTCTGTGCGCCAGGCTACAGTGCAATGGCGTGATCTTGGCTCACTGCAACCTCC

OTI All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence [imiRBase](#). The
Disclaimer: flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms
 should not affect the function of the mature miRNA.

E. coli Kanamycin (25 ug/mL)

Selection:

Mammalian Neomycin

Cell

Selection:

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Both mature and minor (star) miRNA can be over expressed from the construct.

Contains:

RefSeq: [MI0006414](#)



[View online »](#)