

Product datasheet for **SC400545**

MIR567 Human MicroRNA Expression Plasmid (MI0003573)

Product data:

Product Type: MicroRNA Expression Plasmids

Species: Human

Vector: pCMV-MIR (PCMVMIR)

Precursor: GGAAUUCUUAUAGGACAGUAUGUUCUCCAGGACAGAACAUUCUUUGCUAUUUUGUACUGGAAGAACAUGCAAACUAAAAAAAAAAAAAGUUUUGCU

Mature Sequence: AGUAUGUUCUCCAGGACAGAAC

Full insert: >MI0003573 ATATCCCTGGCACACCTTTGTGGGGCAATGCTGAGAACATGAAGAAGGAATGTGTGGAGAAAAGGGG
AGGCTTCCCTGGCTTGCCTTAAGGGATGAAGGCTGGGGTTGGGACAGATGGTGGGTC AAGCCTGGGATCT
TCTCAGCATCTTGACTGGAGATGGGCACCTCTACCTCTACCTCTGATGGCCTCTCCAGTGGATTCCCTAG
AGGATAGGTACTAGGTGGAAGGAGGAGGTTAGGATCACCTCAACTTTTTGTGATGAATCTATGGAAGAGG
ATTCTTATAGGACAGTATGTTCTTCCAGGACAGAACATTCTTTGCTATTTTGTACTGGAAGAACATGCAA
AACTAAAAAAAAAAAAAGTTATTGCTTAAATTCATTTCTTGAGTTGCAGTATGATGCTGTTTATGTTTAA
ATGTCTCCATTGCCTTTTCAGTGGTGAAAATGCCACAGTAACGTATGACCTGCAATTTGGGGTTCCATCA
GATGATGAAAATTTATGAAGTATATGATGAGTGAGGAGTTGGTGTGGGCATTTTGCTACAGGATTTCC
GTGATCAGAATATACCTGGTTGTGAGAGTCTGGGGCTTGATCCAACATCCCTCTTGCTCTATGGTAAGTA
GAGCAAGTAAAGAAGTCAGAGTTGTGACAATAAACAGGTGAAC

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](http://www.mirbase.org). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0003573](https://www.ncbi.nlm.nih.gov/nuccore/MI0003573)



[View online »](#)