

Product datasheet for **SC400564**

MIR586 Human MicroRNA Expression Plasmid (MI0003594)

Product data:

Product Type: MicroRNA Expression Plasmids

Species: Human

Vector: pCMV-MIR (PCMVMIR)

Precursor: AUGGGGUAAAACCAUUAUGCAUUGUUAUUUUUAGGUCCCAUACAUGUGGGCCUAAAAUACAAUGCAUAAUGGUUUUUCACUCUUUAUCUUCUUAU

Mature Sequence: UAUGCAUUGUUAUUUUUAGGUCC

Full insert: >MI0003594 CAGACAAAAGATAAGTTGAGAACTTTGAGAGATTTTTTATTCTAGAAATGCAGTGATGGAACATTAT
CAAGTTTGGGTATTAAACAGTATTCAGAGTAAGAAGAAATCTGAATGCTTTTGTAAAGTTTACAGTCTGA
ATAGATGCTAATTGCAGATGAGAAATACAGCTTTAAGTTTCTTACGTTGTATGTAATTTGGGTTTTTTTG
GTGATTAATATGGCAGTTAGACTCACTTTACACTTTGTAAATTCAAAGAAAAGCAATTATCTAACAAAGT
CAAATTCTTATTAAGACATGGGGTAAAACATTATGCATTGTATTTTAGGTCCCAATACATGTGGGC
CCTAAAAATACAATGCATAATGGTTTTCTACTCTTATCTTCTTATTTATGTAAATGAAAAATAGTAT
TGTCATTTTTTCTGCAATACCACAAAACCTTTGTCTTTCTCATAGTTTACCTTAATTTGTTACATT
GCATATTTCAATTTGAATTCTTATTTTTTAATTAGTTTTAAAAAGCATGTATTATAGAGGATATATGAAG
TATTCATCCAACCAATAGCATAATTTAACTGATCTTTGAAAAATGTTTATTGGTTAAATGCACTTAGTA
ACAATGTATGTAAAGAAGAAACACCATGTAAAACAGGAGCCACTTA

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0003594](#)



[View online »](#)