

Product datasheet for **SC400567**

MIR589 Human MicroRNA Expression Plasmid (MI0003599)

Product data:

Product	MicroRNA Expression Plasmids
Type:	
Species:	Human
Vector:	pCMV-MIR (PCMVMIR)
Precursor:	UCCAGCCUGUGCCCAGCAGCCCCUGAGAACCACGUCUGCUCUGAGCUGGGUACUGCCUGUUCAGAACAAUAGCCGGUCCAGACGCUGCCAGCUGGCC
Mature Sequence:	UGAGAACCACGUCUGCUCUGAG
Full insert:	>MI0003599 CGACACAGTTGTACCCACTTAATGGATGAAGAAAGTAAGGCCAGAGAGTGCATGAACTTGCCTAGGTCA CATGGCTGCTTGTCTCCGGCTCCCGAGGTCAGCGTGATTTTCATGGTCCCGATGGCTCCGGAGTCTGG GACGGAATCCCGTGGCCAAGGGCTTCTCAGGAGCAGGGCCATGGAGAAGGCGGGAGAGGCCTGCGG GTTCTGTGGGTTTTGTGGGTTTTGTGGGCCAGCCTGAGAGACCGACCCTGCTCATAATGCCATCTCTGTG GACAACTCCAGCCTGTGCCAGCAGCCCCTGAGAACCACGTCTGCTCTGAGCTGGGTACTGCCTGTTCA GAACAAATGCCGGTTCCAGACGCTGCCAGCTGGCCCCGTTCTCTGGATTCTGCCTTCTGCGCAAGG GAGCCCGGACGCTGCACTGGAAGGAGCAGGTGGCACTCGCTTGCCCGTGGGCAGGAGCTGTGCAGAGGAA ACGGCGCTGGGATCCCTGGAGGCTGGCAGCATGCTCCGAGGCCGCATCCTGGGGCGGCCGGAGTCTCTTT GGTGGCAGAGCTTGAATCCCACGCCAGCCACCTGCTGCCACTGGCAGCCCCACACCGTGACCCCCAGC TTTGTGCTTGTGGGCTCCTGGCAGGGAGGGGTGGTGCAGGAAGGAGTCGCTC
OTI	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in
Disclaimer:	miRBase . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0003599


[View online »](#)