

Product datasheet for **SC400600**

MIR623 Human MicroRNA Expression Plasmid (MI0003637)

Product data:

Product Type: MicroRNA Expression Plasmids

Species: Human

Vector: pCMV-MIR (PCMVMIR)

Precursor: GUACACAGUAGAAGCAUCCUUGCAGGGGUGUUGGUUGCAUCCUAAGCUGUGCUGGAGCUUCCGAUGUACUCUGUAGAUGUCUUUGCACCUUCUG

Mature Sequence: AUCCCUUGCAGGGGUGUUGGU

Full insert: >MI0003637 CATGGCTGAGCCATGTGCATCCCATGGCTCTCTCCAGAGTCCCTTCCTGCACCTCCAGGGCCTCCCCAG
ATGGAAGTGGCTGGGAGCCGTGCTTTGTTGAAACAAAGGGAAAACCTGTCTCTCTTGGGTTTGCAGCTGC
TGTTGAAACCAACTGAACCTGTTTGTAAATTGTACCTATGTCACAGTAAGATCAGCTGTGCTAGGTATG
TCAGGCCTCTGAAAAGAAATGTCATTAAAGTAAACAAGACTTGGGACACATGTCTGTCAGACTGTACACA
GTAGAAGCATCCCTTGACAGGGCTGTTGGGTTGCATCCTAAGCTGTGCTGGAGCTCCCGATGTAAGTCTG
TAGATGTCTTTGCACCTTCTGTCTCATTGCCATCCCCAATCAGTCCCTGGAAACACAGCTGTGGCCCTG
ACCCTTTTAGCCCCCTGCATGGAATGCACAGGCAAGGCGTGTGGCCTTTGGGACTCCAGCCGACTAGAGCA
AGACAGAGGAAGAAGAGACGGGGCAGAGGGCAGCCGATAGTTGGAACTTCTGTCCCTAAGCCAAAAGG
TTCTCAAACTGAGCAACAGGATGATTGTGTCTCCTGATCGTATTGCCATGGAACTGCCTTCTGTTCT
GCCTAGAAGAAATAGTACTGTTTGAAGTATCTCATCCATC

OTI All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in
Disclaimer: [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0003637](#)



[View online »](#)