

Product datasheet for **SC400637**

MIR660 Human MicroRNA Expression Plasmid (MI0003684)

Product data:

Product Type: MicroRNA Expression Plasmids

Species: Human

Vector: pCMV-MIR (PCMVIR)

Precursor: CUGCUCCUUCUCCCAUACCCAUUGCAUAUCGGAGUUGUGAAUUCUAAAACACCUCUGUGCAUGGAUUACAGGAGGGUGAGCCUUGUCAUCGUG

Mature Sequence: UACCCAUUGCAUAUCGGAGUUG

Full insert: >MI0003684 ACGTTCCCAATGCCGAGTCCATGCAGACTGTGGCAGGGAAGCCATCATTTCACTACTGTCTCCACCTCCT
CCATAACCACACATTTCAGTGAGCACATGAGGCTAAGTGAACACACAACAAGTACTTTATTAGGGAATAGA
CTTGCCACTGCTTCTCCAGGCGTGACCTGACACTTTAATCTCAATCTGTCTCTCTCAATCAATCAAT
CAATCTCTCTCCTCTCTCTCTTTCTCTCTTTCTGTCTTTCTTTCTCTCTCCTCCCCACCCCCCTCCC
ATGTGCAGAAGAATGAACTGCTCCTTCTCCCATACCCATTGCATATCGGAGTTGTGAATTCTCAAAACAC
CTCCTGTGTGCATGGATTACAGGAGGGTGAGCCTTGTCATCGTGAATCTTCTGCTGAACAATCCAAGTAG
CCCGTGTTTATTGCGCACTTTTGATGTGCCATGAGCTGTGCTGGGAGTACCCTGCTAAATAATGGACAGT
GTCTCTCTCCGGGACTTCCTATCTAGTGGAGGTATTAGGTGCCTAGACTTCCCCAGGCCTAGGCAACAGA
CACAGGTGTGTGAGGAATTGAGCCTTCTCATTCTTATGAATACACATGCTCACATTGTGCATGCATGT
GCACACACATGTGTAAGTGTAGAAACATTTCTGTACAAAGCTGAAAGGAAGGAGCC

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0003684](#)



[View online »](#)