

Product datasheet for **SC400649**

MIR7-2 Human MicroRNA Expression Plasmid (MI0000264)

Product data:

Product MicroRNA Expression Plasmids

Type:

Species: Human

Vector: pCMV-MIR (PCMVMIR)

Precursor: CUGGAUACAGAGUGGACCGGCGGCCCAUCUGGAAGACUAGUGAUUUUGUUGUCUACUGCGCUCAACAACAAAUCCAGUCUACCUAAUGGUGCCAGCCAA

Mature Sequence: UGGAAGACUAGUGAUUUUGUUGU

Full insert: >MI0000264 CTGGAGTGAGCTGTGATTGCACCACTGCACTACAGCCTGGGTGACAGAGCAAGACCCTGTCTACCAAAA
CAAAACAAAACAACTGTGTCCATTATAGAAAGTAGAACACTCCAGAAAAAGGGAAAAGTAAGAAAAAAT
AATTTCTCAGGGAATTCTGAGGTGAGGAGAGGTGTTCAAAGATGAACCTTAGTAATAGGAGCGGAGACAA
CGCCTCGCCTGAAGGAGCATCCAGACCGCTGACCTGGTGGCGAGGGGAGGGGGTGGTCTCGAACGCCT
TGCAGAACTGGCCTGGATACAGAGTGGACCGGCTGGCCCATCTGGAAGACTAGTGATTTTGTGTTGTC
TACTGCGCTCAACAACAAATCCAGTCTACCTAATGGTGCCAGCCATCGCAGCGGGTGCAGGAAATGG
GGGCAGCCCCCTTTTTGGCTATCCTTCCACGTGTTCTTTTTGTATCTTTTGTGTTTCTAGAAAACAT
CTCAGGTCACCACCCTAAACACCCATCGGCAGACTCTTCCAGAAACATGGACAAACTATGGTGTGAGCAG
GGAAGTGGGGCAGGGTCTGTGTTTTAGGGAGCAGTGGACTGCTTCCATACCACTACCTGAAGGGCAGGG
GGAGAAAAACCTCCTTAATTACTAATTCCCCAACAAATCACTTTGGGGCACCTATAATGGAGGCACCT

OTI All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence [imiRBase](#). The
Disclaimer: flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms.
This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0000264](#)

Synonyms: MIRN7-2 mir-7-2 hsa-mir-7-2

Locus ID: 407044



[View online »](#)