

## Product datasheet for **SC400690**

### MIR936 Human MicroRNA Expression Plasmid (MI0005758)

#### Product data:

**Product Type:** MicroRNA Expression Plasmids

**Species:** Human

**Vector:** pCMV-MIR (PCMVMIR)

**Precursor:** UCAAGGCCACUGGGACAGUAGAGGGAGGAAUCGCAGAAAUCACUCCAGGAGCAACUGAGAGACCUUGCUUCUACUUUACCAGGUCCUGCUGGCCCCAGA

**Mature Sequence:** ACAGUAGAGGGAGGAAUCGCAG

**Full insert:** >MI0005758 GAGCAAGATATGGTCCCAGAAGGCCACCATTTTGGAGAGAGGTTACAATAATGCCTAGAGCAGAGGAGAC  
AAGAACATGGAGGACCTTGACTTCTAGGCCCTCAAGTAAACGCGACAAGCATTGAGCATCGAGGCGCTTG  
AGATGCAATGGGAAGGCACTGAACAATTCTCAACAGGAAAAGAGACAATTGCTGGTTTAGGAGAAGGGGC  
AGTGGGAGTACCGTATCCAGAGAGTGTGACAGACAGGGCTGTCCCGTGGCTGATGGCTGCGTGGGCTCCT  
GGGACACAATGTGCTTCAAGGCCACTGGGACAGTAGAGGGAGGAATCGCAGAAATCACTCCAGGAGCAAC  
TGAGAGACCTTGCTTCTACTTTACCAGGTCTGCTGGCCCAGACGGACACCAAGGCCAAGAGGTTGGTC  
ACTCATCTGCTGCCCCAGCCTCACTCCCGTCACAGCCCCTGTTCCAGCATGCAGCCCTCAGTCCCCTTA  
ACCATAGTGCAGAGGAGCATATAGCATGTGGAAGTGTGCAGTGTCCAGTTCGTTCCCTTCTCCTTTCC  
ACCTGCCAGCAAGTCAATTAGGAAAACCCTCAGAACTATCAGTTCAGATGTGATAACCAAGTTTAGCTTT  
CCTAAATGTAAACACTTCCACGTTGGCCATTCTCTCTAGGTGAACAAGGTTACTGGG

**OTI Disclaimer:** All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

**E. coli Selection:** Kanamycin (25 ug/mL)

**Mammalian Cell Selection:** Neomycin

**Components:** 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

**Also Contains:** Both mature and minor (star) miRNA can be over expressed from the construct.

**RefSeq:** [MI0005758](#)



[View online »](#)