

Product datasheet for **SC400703**

MIR96 Human MicroRNA Expression Plasmid (MI0000098)

Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Human
Vector:	pCMV-MIR (PCMV MIR)
Precursor:	UGGCCGAUUUUJGGCACUAGCACAUUUUUJGCUUGUGUCUCUCCGCUCUGAGCAAUCAUGUGCAGUGCCAAUAUGGGAAA
Mature Sequence:	UUUGGCACUAGCACAUUUUUJGCU
Full insert:	>MI0000098 TAGGAGCAGGGAACGGGCATCGTGGGCCGCTGGTCTCTCCGAGGGTCGGCAGGCCGAGAGTGTGACTC CTGTTCTGTGATGGCACTGGTAGAATCACTGTGAACAGTCTCAGTCAGTGAATTACCGAAGGGCCATA AACAGAGCAGAGACAGATCCACGAGGGCCTCCGGAGCACCTTACCCACTTCTGCCTTGAGTGCTCCTAGA CGTCGGAAACAGGCTGCTTCCAAGGGTGCAGGGATGCAAGGCCCTCGTCCAGTGTGTCAGAGAGCC CGCACCAGTGCCATCTGCTTGGCCGATTTTGGCACTAGCACATTTTGTGTGTCTCTCCGCTCTGAGC AATCATGTGCAGTGCCAAATATGGGAAAAGCAGGACCCGAGCTGCGTCCGCTCCCTGCATCCTTGTGT CAGGGCCCCAGCCTGCTCCTCAAGGCCTCCTACCGCCTCCCCAGCCATCTGGCTCAGTGTGTG TGAGGGCCCAGCGCTGGTGGGCAGCCAGATCGCCTTACACTGCCTGGGGCCACGGTAGAGCTGGGAGCCC AGCAATCTGAGCTGGGCCAGCAGATGGGGCCGCCAGGGCAGAGGTGGGGGAGTCTGAAACCATCTGTAG GCCATCCTGAATGGTGCCGTGGTTGGAAA
OTI Disclaimer:	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in miRBase . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	MI0000273 Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0000098
Synonyms:	DFNA50 MIRN96 miR-96 miRNA96 hsa-mir-96
Locus ID:	407053



[View online »](#)