

Product datasheet for **SC400722**

Mir103-2 Mouse MicroRNA Expression Plasmid (MI0000588)

Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Mouse
Vector:	pCMV-MIR (PCMVMIR)
Precursor:	GUCUUCGUGCUUUCAGCUUUUACAGUGCUGCCUUGUAGCAUUCAGGUCAAGCAGCAUUGUACAGGGCUAUGAAAGAACCAAGAA
Mature Sequence:	AGCAGCAUUGUACAGGGCUAUGA
Full insert:	>MI0000588 CCTCCTCTTCATTGATTGTCCTCCCATTCACCATACTCCTTGCTCTTCTGAAGATAAGCACTACT GCAAGAGGTGAATGGTTTCTTTGAGGAGACTCTAGAGGGGGCCAGGAAGAAGGCCAGTCTTGAGCTACTG ATGGGATGAATGACGGGTCAATGGTGACTAGACGCTGTGTTTGCTGGGTGCTGTTGAGGTGTGGACAG GAAGCCACAGGAAGCCACTGGAGCTGGGAGAAGACTGAGCGGTGAGAGGAAGGCCACCTTGGTTCTTCT CCACCTCAGCAGAGCTGTCTTCGTGCTTTCAGCTTCTTTACAGTGCTGCCTTGTAGCATTGAGTCAAGC AGCATTGTACAGGGCTATGAAAGAACCAAGAATGGGCTGCCCTTAGGGAAAGAGGATCTAGGATCCAGTT GGTGCAGCTTATGGCTTGTCCCTATTTATACACTTTGGTTTACAATAATTCTTTAATAATATAGGCATG TGCCTTAAAAAGAGACCTCAAGACCTAAAAGGATTTGTCTAAATTTAATTTATTTAGAAGTTATTAAC ACCAGATTTTAAAAATCCATCTTCATTTGAATAACATCACTTTATATGTTTGCTTTCCTGGGCTTAAG TATCAAAAACCTTTTAGTTGTTGTGCATACATTTTATGC
OTI Disclaimer:	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in miRBase . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0000588
Synonyms:	Mirn103-2 mir-103-2
Locus ID:	723825



[View online »](#)