

Product datasheet for **SC400732**

Mir1188 Mouse MicroRNA Expression Plasmid (MI0006290)

Product data:

Product MicroRNA Expression Plasmids

Type:

Species: Mouse

Vector: pCMV-MIR (PCMVMIR)

Precursor: AUACUCACAGUCUCCCAGCUGGUGUGAGGUUGGGCCAGGAUGAAACCAAGGCUCUCCGAGGCUCCCCACCACACCCUGCUGCUGAAGACUGCCUAGCAAGGCUG

Mature Sequence: UGGUGUGAGGUUGGGCCAGGA

Full insert: >MI0006290 ATCGCTCGGTCTGTCAGTCAGTCGGTCGGTCGATCGGTTCGGTTCAGTCGGCTTCCTGTCTTCGAACT
GCAATCATTCTTATGAGAAAACAAACAGGAAATACTTTGCCACACACCAATGCAATCAGGGTGGGTGATG
GCAGTCGGGATTCAAGTCACCCTAGTTGAAGCCCTGAGCTGGGTGGGCTGTCTGTCCAGAACTCACATGG
GAAATTTACAAAGAGGAAAAGTCTCATGATCTCCATTGGATTGGTGAAGACACCATGGCACAGAAAAGC
AAAGAAGTCCAATACTCACAGTCTCCAGCTGGTGTGAGGTTGGGCCAGGATGAAACCAAGGCTCTCCG
AGGCTCCCCACCACACCCTGCTGCTGAAGACTGCCTAGCAAGGCTGTGCCGAGTGGTGTGGGTGGGCAAA
GGATGGATTGACTAGTTATAAGCACCTCCCAGAAGTCCTTGAGATCTCAGCTATTTTGTGGCCAATCCC
TTTGAAGGGCAGGGTCAAAATTAGGTGTAGGCACCTACCTATTAGCGGAGGCTTATGAAAATGACATTAT
TAACCTGATGCAGAGAGGGCCCTATTCATCTGTAATTAATTAATAATAGAGTTAGTGGGGGGCTAGGGT
CTCTGTGATTGACTGCAGTCTAGAGTGCTCAGAATTTCCAGGGCTCTGGTGAACCTTTTC

OTI All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence [imiRBase](#). The
Disclaimer: sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect
the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0006290](#)



[View online »](#)