

Product datasheet for **SC400735**

Mir1192 Mouse MicroRNA Expression Plasmid (MI0006297)

Product data:

Product	MicroRNA Expression Plasmids
Type:	
Species:	Mouse
Vector:	pCMV-MIR (PCMVMIR)
Precursor:	CCUCCAAACAGUAAAAACAAACAAACAGACCAAAUUUGUCAUUGUUGUUCGAUUUGCUCUUUUUGUUUGUUUGUUCUUUGAGACAGGGUCUGG
Mature Sequence:	AAACAAACAAACAGACCAAAUU
Full insert:	>MI0006297 AACACAGGAGGAATCAGACTCGTCCAAGGTTATGCTGAGATATTCCTACTGTTGTTCAAGGTCTAAGTTG TACGTGAGGGCTGTTGAGTGGCTTGCTCCCAGTGAGAAGAGATTCTGTATTTATTCTTATGGCTCTTTTG GGTTTGGCCCTAGAGCCACCATCTTAAGCTGTCTTCCTTGTCGCCAGGTCTGATGTGAAGATAGTTTTT ATGAGTTTTAGTTAGGAAAAGTCAACAAAATAAAGTCTCCTGACCTTGTACCCGAACTATTTCTTCC AAACAGTAAAAACAAACAAACAAACAGACCAAAATTTGTCATTGTTGTTTCGATTGCTCTTTTGTGTTT TTGCTTTTGTGTTTGTAGACAGGGTCTGGCTTGATAGCTCATGCCAGCCTGCAACTCTATGGTCTAGGC TGCTTTTCTGCCTCAGCTTCCAGACTGCTGAGATGACAGGCCAGTGCTGCCACACCTGCCAGGGCTGGG TTTGGGATCCATAAATAGGTTGCAAGCAAGAAAGTCCAGACTTGACTTGCTCTGAGCTGGGCGTTGTGT CTCTCACCTATAATCCCAACACTTGATAGGCTTAGGCAGCAGACAACTGCCACGTGTTAGAGGGCACACT AGGCTACTTAGCAAACTTAAAAAAGATTAAAAATGGAGATCCAAGTTCAGGAGTTGC
OTI	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence imiRBase . The
Disclaimer:	sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	<u>MI0006297</u>



[View online »](#)