

Product datasheet for **SC400747**

Mir124-2 Mouse MicroRNA Expression Plasmid (MI000717)

Product data:

| | |
|----------------------------------|--|
| Product | MicroRNA Expression Plasmids |
| Type: | |
| Species: | Mouse |
| Vector: | pCMV-MIR (PCMVMIR) |
| Precursor: | AUCAAGAUCAAGACUCUGCUCUCCGUGUUCACAGCGGACCUUGAUUUAAUGUCAUACAAUUAAGGCACGCGGUGAAUGCCAAGAGCGGAGCCUACGGCUGCACU |
| Mature Sequence: | CGUGUUCACAGCGGACCUUGAU |
| Full insert: | >MI000717 CTTTCGATCCATGTCATCATTTGAAAAGTACTAGATTTATAGACATATGTATGGGTTTTAGGTGCGCTGTAAT GGCATGCAGATGTATGCATAAGAATCGCAGACTCATGCACCCGCTACTCTTCAACCTTGCAGACTAATTA ATGGCAACTCTGACAGATGTGGCGGGAGGAAAGCTGAGGTGGGAGTACTGCTCAGAGCTACAACCTTAGG AGTAGGGACTCCAAGCCTAGAGCTCCAAGAGAGGGTGAAGGGCAGGGAGAAAATTATAGTAATAGTTGCA ATGAGTCACTTGCTTCTAGATCAAGATCAGAGACTCTGCTCTCCGTGTTACAGCGGACCTTGATTTAAT GTCATACAATTAAGGCACGCGTGAATGCCAAGAGCGGAGCCTACGGCTGCACTTGAAGGACATCCGAGA GAAGTTAGGAAGGGTGGGGAGAAACAATTCTAGAATGAACCCATCCTGTGCGACACAAGGACCCCTCGCA TTGATTCAAAATCTTTAGAATACAAGACAGTCAAATGTTTTTCAGCGAGAATCTGGGAAATCCGCATGCA TACTCCCAAGGGGAATATTAGAGGTTTCATTGACAGCTGTCAGGAGAGCTGTAGTTCCCATTTGGTGCAA CGGTCTAAGAGATGGCACTCAGAGAGTGTGGCCTTTTCAAATAAAATCTGGGTGCAACTGGGAAAAC AA |
| OTI | All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence imiRBase . The |
| Disclaimer: | flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA. |
| E. coli Selection: | Kanamycin (25 ug/mL) |
| Mammalian Cell Selection: | Neomycin |
| Components: | 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube |
| Also Contains: | Both mature and minor (star) miRNA can be over expressed from the construct. |
| RefSeq: | MI000717 |
| Synonyms: | Mirn124a-2 mir-124a-2 |



[View online »](#)

Locus ID: 723950