

Product datasheet for **SC400748**

Mir124-3 Mouse MicroRNA Expression Plasmid (MI0000150)

Product data:

Product Type: MicroRNA Expression Plasmids
Species: Mouse
Vector: pCMV-MIR (PCMVIR)
Precursor: CUCUGCGUJUCACAGCGGACCUUGAUUUUUAUGUCUAUACAAUUAAGGCACGCGGUGAAUGCCAAGAG
Mature Sequence: CGUGUUCACAGCGGACCUUGAU
Full insert: >MI0000150 CGGCTCTGCATCTCCCTCGCCGCGGCTCCGAGGGTCCCTACTGCAGGAACGCCCGAGGGATTGCCCGGGCCCTTCTCTACCTCGGGGAGGGGAAGGAAACTGGCTCCTCGGAGTTGGGGGCGGACCAGGGCAGCACTGCCAGCTTAGCGCGGAGCCCCGCCCGCCCGCCCGCCCGCGGCGAAGAAGCTGGAGCATTGCGGCCCTCAGGCTAGGACCCACGCAAGCCGCAGCCCGGCCCGACCCGGTCCAGCCGCTCGCCGCCCGCCGAGCCCTGAGGGCCCTCTGCGTGTTCACAGCGGACCTTGATTTAATGTCTATACAATTAAGGCACGCGGTGATGCCAAGAGAGGCGCCTCCGCCGCTCCTTTCTCATGGAATGGCCCGGAGCCCGCCGTCGCGGCGCCCTCCCGCGGGAAAAAGGCGAGCCCGGCCCGGTCGCTCCGCGGACAAATCCGGCGAACAATGCGCCCGCCAGAGGGCGCCAGCTGCTGGCCCGGGATCTGGCCCGGGACACAAAGGGGCCCGCATGCCTAGGACGTACGAGAAGGGTGGGGGCGACCGCGGGGAGAATGTGGCCACCCGGCCCTGACCACGTGTGTCGCGAGCTCTGAGTGCT

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0000150](#)

Synonyms: Mirn124a-3 mir-124a-3

Locus ID: 723951



[View online »](#)