

Product datasheet for **SC400793**

Mir155 Mouse MicroRNA Expression Plasmid (MI0000177)

Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Mouse
Vector:	pCMV-MIR (PCMV MIR)
Precursor:	CUGUUA AUGCUAAUUGUGAUAGGGUUUUGGCCUCUGACUGACUCCUACCUUGUAGCAUUAACAG
Mature Sequence:	UUA AUGCUAAUUGUGAUAGGGU
Full insert:	>MI0000177 AATATTTCCACAGAATTTAAAGTTTGATTTTATAAATTAGTTTTAAATTGCAAACCTTACAAAAATTTTC CAAAATTTACTTAGAAGCATTTTCAGAGCTCTTTTCTTTCAAAGCTGAAGTCTACCTTGCCTTCAGATTCA TAATCGTTCACCTACTTTATAACTAAGTGTGTGCGTACACACAACTATTTCTTTTCTCTTTGCAGGTG CTGCAAACCAGGAAGGGGAAGTGTGGTTTAAAGTTGCATATCCCTTATCCTCTGGCTGCTGGAGGCTTG CTGAAGGCTGTATGCTGTTAATGCTAATTGTGATAGGGGTTTTGGCCTCTGACTGACTCCTACCTGTTAG CATTAACAGGACACAAGGCTGTACTAGCACTCACATGGAACAAATGGCCACCGTGGGAGGATGACAAG TCCAAGAGTCACCCTGCTGGATGAACGTAGATGTCAGACTCTATCATTTAATGTGCTAGTCATAACCTGG TACTAGGATAGTCCACTGTAAGTGTTACGATAAATGTCATTTAAAAGATAGATCAGCAGTATCCTAAAC AACATCTCAACTCAAGCCACATGTTTATTTTTATCTTGAATGGAAAGTAAAAGTGTATCATTTTTTA TTTCAAATATGTTTCATAACCATCTT
OTI Disclaimer:	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in miRBase . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0000177
Synonyms:	Mir155 mmu-mir-155
Locus ID:	387173



[View online »](#)