

Product datasheet for **SC400809**

Mir185 Mouse MicroRNA Expression Plasmid (MI0000227)

Product data:

Product Type: MicroRNA Expression Plasmids
Species: Mouse
Vector: pCMV-MIR (PCMVIR)
Precursor: AGGGAUUGGAGAGAAAGGCAGUUCUGAUGGUCCCCUCCAGGGGUGGCUUUCUCUGGUCCUU
Mature Sequence: UGGAGAGAAAGGCAGUUCUGA
Full insert: >MI0000227 TTTGCTCTCAGCCACAGGTGAACAGGAGCAAGAAAGGGGCTAAAGAAGTGGCATAGCAGAAGGGACCTC
ACCAAGTCCCGAAAAGAGCATCAGACCCTCGGGTGCCGTTGGAGATGCCTGGCATTACAGGGGCTTGC
CAGCCGCTTCTGGGAAAACATCTGGGCTTAGGTAAAGGCAGGATCTGTTACGGGAGAGGTATATCAGGA
AAAAGGTATGAGGAGTAATTTGCTGGCCAGAGACGGCAGAGGGACTGTGGGCAGCGGGCTGGCTTGAG
TGGGGGTGAGGGATTGGAGAGAAAGGCAGTTCCTGATGGTCCCCTCCAGGGGCTGGCTTTCCTCTGGT
CCTTCTCTCCAATGACTGCGTCTTCATCGATGCTGCAACCCTGAGCTCTGGCAGCACCTCCAACCTCCAG
GCTGCCTGGAGACACATCGATGGCCTTCTGTGGACCTTGGAGCCCTCCTCACATCTCCAACCTAGAGGG
CTGGTCTGCTAGGGCCAGCTAGGATGGCAGTCTGCAGCTGAACTAGAGAAAACCCCTTGATCTGGCTA
TGTCTGCTGATAGCTGCCTCTGCCAGAGGTCTACCCAGTACAGCCTAGGGAGGACCAAAGGAGTCTAAG
AGCCA

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)
Mammalian Cell Selection: Neomycin
Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq: [MI0000227](#)
Synonyms: Mirn185 mmu-mir-185
Locus ID: 387180



[View online »](#)