

## Product datasheet for **SC400810**

### Mir186 Mouse MicroRNA Expression Plasmid (MI0000228)

#### Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Mouse
Vector:	pCMV-MIR (PCMVIR)
Precursor:	ACUUUCCAAAGAAUUCUCCUUUUGGGCUUUCUCAUUUUUUUUUUAAGCCCUAAGGUGAAUUUUUUGGGAAGU
Mature Sequence:	CAAAGAAUUCUCCUUUUGGGCU
Full insert:	>MI0000228 CGTTTCTAACAACCTTCTGTAGAAAGACTCACTTACTTGAGCCATTCTAACATTTGAAGGCAGCTGC CATTGAGTCCCCCTGAGGAGAGAGAGGGTAATGGAGCACTGACCACTGCTCTGCTGACCAAGCGCTGAGT AGAGCAAAGGTTTCTTTCTTATTGCATTAACATCCGTGCTCACTTCAGACTGTGGAGTCTTGGTGTGGC CCTCTGGGAAGTGTAGGACTGAGGTTTACAGGACAGCCTTCATATCCTTCCAAACATTTAATTGCTTAC AACTTTCCAAAGAATTCTCCTTTTGGGCTTTCTCATTTATTTAAGCCCTAAGGTGAATTTTTTGGGAA GTTTGAGCTAAATGCCTTCAACTAAAAATATGTAATAGAGATAAATTTGCTCTTTAAATTACACATTTA CTTCCACTTAAAAATGTAATGTTACTACTTCATACTTCATAGCACAGTTCTCATACAAATCACATTGT AGAGTTGCCATTTTCATTGACTACTTCACAGCTTCAGTGTCTTTTTCTGGGCCATGTCCTCAGGTGTC TGTAGCACCTGTAAGAAGTTAAGTTCATTGTCTTCAGCAAGCTTTATAAAAAGTGAATTTCTGACAAA TGGGAGG
OTI Disclaimer:	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in <a href="#">miRBase</a> . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	<a href="#">MI0000228</a>
Synonyms:	Mirn186 mmu-mir-186
Locus ID:	387181



[View online »](#)