

Product datasheet for **SC400813**

Mir1892 Mouse MicroRNA Expression Plasmid (MI0008320)

Product data:

Product Type: MicroRNA Expression Plasmids
Species: Mouse
Vector: pCMV-MIR (PCMVMIR)
Precursor: UGGAUUAUUUUCUACUCUUGACCAAAGUUUUGUAAGAAACAUAACAGAAUUUGGGACGGGAGGGAGGAUACUUUGAGAA
Mature Sequence: AUUUGGGACGGGAGGGAGGAU

Full insert: >MI0008320 GAGTAGAACAGGAAAAATCATCCTGGCTCATGTGTTGTGTTCTTTCTTTTATTTTAAACAGAGGCTCTT
 ATATAGTAGCTTTTATCTATTTTAACATTGTAGTCATTTGTACTTTGATATCAGTATTTTCTTAACCTTT
 GTGACTGTTTCAATATTATCCAGTGAAAGCTTTTCTTAATGTAACCTTTGAGTACATCTCAATTGCCTTCT
 ATTTTAAAAACCTAAGGTCATTAGTTGGGCTTTACTGTTCTTGCTATCATATGGCATATACATCTGCCTG
 GATATATTTCTACTCTTGACCAAAGTTTTGTAAGAAACAATACAGAAATTTGGGGACGGGAGGGAGGATAC
 TTTGAGAACTACTTACAAAAGATTTATCTGCTTGCTTGAACCTCAGGAGTACGGTTTTAGCTCCCTAGACT
 CTTAACAGCTTTTGCTCTAAACTATTAAAGTGTTCCTTAGTAATGAAAAAGTAAAGATCTTGCTAACAT
 TAATACAAGGAACATTTACCTTTTAGATATTAAATCTTATGTGGACTTATTTCCAACCACTTTGATG
 ATTCCTGTGATTGGTGATTGGTGGCCGGCTTGATGACATATGCTGAGATACACTGTGGAGCTTTTCATA
 GAAGGCTTATCAGATCCTCTGAGG

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0008320](#)


[View online »](#)