

Product datasheet for **SC400829**

Mir1905 Mouse MicroRNA Expression Plasmid (MI0008315)

Product data:

Product Type: MicroRNA Expression Plasmids
Species: Mouse
Vector: pCMV-MIR (PCMVIR)
Precursor: GUGUUGUCUGCUGCUGGAUGCGUUUGAUGGUCGCUCCUUGGGGCCACCACCAGUCCACCACGCGGUAGGGUACUCGUA
Mature Sequence: CACCAGUCCACCACGCGGUAG

Full insert: >MI0008315 AGTCCACGCCGATTGCCTATGCAGCCCAGGCTGTTCTGCCTGAAGGTGGAGAGCGGCTTACAGCCTGG
 AGCATGCACCCGCCAGGCGTCTGAGTAGCGGCTGTCCAAGGCAGCGTCAGGGCTGCCAGCCAGGAAGTCG
 CCGTCGCTGTTGTATTCAAGAATCTTGCCGGTGCGCACAGCGATGTGCGTCTCGATCTCCTCCCGCGCAC
 GCTCCACGTTGCCCGGTGCGCCGTTATCTCAACACGGGGTCGCGGTCCCGGCTGGGAGTGATGATGTA
 CGTGTTGGTCTGCTGCTGGATGCGTTTGATGGTCGCTCCCTTGGGGGCCACCACCAGTCCCACCACGCGG
 TAGGGTACTCGTACACGGATGGTCACTTGGCCTGGCAGAGCAGGGGCCACGCCAAAGCGGCGCCAGACT
 TGTTGCGTGAGGCTCTTATCATGGAGAAGTGCTCAGCTGCAGAGATGATCTCCCGCCGGGCTGTGGCCAC
 GTCTTCCCGCCTCCAGTCACCATGAACACTGGCTCCTCACCTCGCACAGGTGTCTTGATGTAGGTGTTG
 GTTTTGGCCCTCAGAGCCTTAATCTTGACGCTGGACAGAGTGGGAAAAATTTACAGGACAAAGGAGT
 TTGTACCTTCCCCAGAAT

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0008315](#)



[View online »](#)