

Product datasheet for **SC400871**

Mir1955 Mouse MicroRNA Expression Plasmid (MI0009950)

Product data:

Product Type: MicroRNA Expression Plasmids

Species: Mouse

Vector: pCMV-MIR (PCMVMIR)

Precursor: CUGGUUUACUUACUACAAGUCCCAGGAUGCACUGCAGCUUUUUCUUAACUUCAGAUAAAGAGCAUUGCAUGCUGGGACAUGUAGUUGGUUUAAAUAGCA

Mature Sequence: AGUCCCAGGAUGCACUGCAGCUUUU

Full insert: >MI0009950 AGTGGGGAGATGGTAAAAATGTTGCAACTGATCCAACAAAACCTATTTTTAAAAGTATAAAATTATGTTT
CTGTCATGTCAGTTATTTTTAAAGCATTTTGTCTTACTTGAGGTTTTCTAGTGAAGTTATTAATGCG
AATTGCTACTTTTTCTAATATGTACTAGAGAAACAATACCAGAGCAGAGAGCAAACAAATTTGATTTGAA
TTTGATTAGACTCTGAAAATAAGACCTAATCAGATGCGCTCCACAGTTCAGACTTTCTGAACTTGCTGT
GAAAAGCTGGTTTACTTACTACAAGTCCCAGGATGCACTGCAGCTTTTCTTAACCTCAGATAAGAGCAT
TGCATGCTGGGACATGTAGTTGGTTAAATAGCACGGAGAACTCTTTTTCTGTCTACTCCATTTTATAA
ATAGATCCCTAAACTCTAAGAACTTCATTTTTTAAAGTTACATTAATGAAATCATACTTACTCTGTTT
CTGAAATTATAAGTGAGACTTTATAACAGCTTGGTGTCTGTTTATCACAGTGCTTTTACTAGCTTT
TATAAGATCTAAATGCAAAATAACATATTCAAACTTGACTTGAACTGTAAGTAGTTTCTGTGCATGT
GGTCGAATAGGTTATTGCAGAAGTCTGGTATCATTGTTTTCTAGTGAGG

OTI All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in
Disclaimer: [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0009950](#)



[View online »](#)