

Product datasheet for **SC400880**

Mir1964 Mouse MicroRNA Expression Plasmid (MI0009961)

Product data:

Product Type: MicroRNA Expression Plasmids
Species: Mouse
Vector: pCMV-MIR (PCMVMIR)
Precursor: CUUUCUUGCGAAGCUGGAGCACAAAAGCCGGUGCCGGUGUGAUGCCUACACCGACUUCUGGGCUCCGGCUUUGGAAGGAAGGA
Mature Sequence: CCGACUUCUGGGCUCCGGCUUU

Full insert: >MI0009961 AACTTTTATATAAAAAACAATCATTTTCCTTTAACTTTTTCTAAGCCCTCAGGGTGCACCCTACACAT
 TAAATTTTTTTACTAAATCATAGGATCAGAAAGCCGAGGGCATAAATGGCTACAAGAAATTAACCATTAC
 CCTTAATCACCCACCCAGCTTGATCAAGAAAGGCACCTCAGCTAGTAATAACCAGGCTTACTCAACTTTC
 AAAAATGAACCTTTCTGAACCTTTAAATTTCTGCCAAAATTTCTACATCAAAGCAACTAACACAAACAT
 CTTAACCTTTCTTGCGAAGCTGGAGCACAAAAGCCGGTGCCGGTGTGATGCCTACACCGACTTCTGGGCT
 CCGGCTTTGGAAGGAAGGAGGTAGCCTCGGCTGGCATTGCTGGTATTGCCTTGACTAATCCTCACACCTA
 CTTGCAGAGGTATGCTAGTACCACAGGGAGGAGGGACAGATGCCAAAGGGAGCAGGGAGTCAACTGTCC
 GTCTTGCTTTTAGTTCTCCTACAGGCAGTGCTTTTTTATTTTCTAAATGCCTCTGTGTCCTTCTCCCC
 AGGAGTATGCCCTTGCTATTTTATGTCTCCCTATTACTATTTTAATCATAACCAAAGTGCTTCCACTTG
 AGCGTGTCAAGTTCGGAAAACAGGTGGC

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0009961](#)


[View online »](#)