

Product datasheet for **SC400883**

Mir1967 Mouse MicroRNA Expression Plasmid (MI0009964)

Product data:

Product Type: MicroRNA Expression Plasmids
Species: Mouse
Vector: pCMV-MIR (PCMVMIR)
Precursor: CCUGAAUUCUGUUUCUUCUCUCUUUGCUCCUUUGUCUGACAUUCAUGAGGAUCCUGGGGAGAAGAUGCAGAAUCCUGG
Mature Sequence: UGAGGAUCCUGGGGAGAAGAUGC

Full insert: >MI0009964 GCCTGTAGTAAGCTTTCCAAAGTGCACAGTCTCTTGCCTCAACACAGCAGCGTGGTGAATGTTCTGTCAG
 GGTTTAAGCTAAGTTGTTTACCATAATGGAAACAACTTTATGACTTTTAAAACTCAAAGGTTTAGCT
 GGGTTTCATAGCAAGTGTAGGTTGTTAATGTATCCTTTCTTTGTTACTTGGGATTTTAAAGATTGTATA
 TGGTTTTTGTCTACTGTAGCTGTGTGTGACGAGTTACTGTTACAGTAGAAGTAAGAGCATTATTTATTGAA
 ATTAATGCCTGAATTCGTTTCTTCTCTTTGCTCCTTTGTCTGACATTCATGAGGATCCTGGGGAG
 AAGATGCAGAATTCCTGGAAATTGACTTACAGATCGTTGTGGGTGCTAGGAAGTGGACTCGGGCTCTCT
 ACAGGAGTCACAAGTGTCTTTACTGCTGAGCCGTTTCCTCAGCCACATCACCTAATTTCTTTTCACTT
 CAAGAAGCATTATTTCTTCAAGTACTGTATTGCTTTTGTCTTATTTAATAGTGCTACTCCTGTT
 ACACAACACAGTTTGAAGCTGTGCTAAACCATCTCTTCGATCTTAAATCATGAATGGTGACTGTCATC
 CTGTGAGTGATTACAC

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0009964](#)


[View online »](#)