

Product datasheet for **SC400884**

Mir1968 Mouse MicroRNA Expression Plasmid (MI0009965)

Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Mouse
Vector:	pCMV-MIR (PCMVIR)
Precursor:	CCUCCUCUGCAGCUGUUAAGGAUGGUGACUUAUUGUUAGGGCCACCACCUCUGCAGCUGUUAAGGAGG
Mature Sequence:	UGCAGCUGUUAAGGAUGGUGACU
Full insert:	>MI0009965 GTTGAGTGGGCACAGAGTAGGAGGAGGTGGGGCCTGAGGCGCGCTGTGGAGGCTTCCTGTGGGTGTGGG GGAGCATGAAGGTACTGACCTCTCTCGTGTTCCTGGACGCCTGCTGTCATGGATCCCACCTGCAGATTTT GGTCTTGTGGATTGCCAGTTGGGGAACAGTTTCCTGACTGTGTTCAGGGTGTGCTGGGGCCTGTGAG CTTGCAAAGTGCATGGCTTAAGTATAATGATGCTGGAATGTGAGACCCCGGTATGGGTGGATCTACTG TCAGGGCCACCTCCTCTGCAGCTGTTAAGGATGGTGGACTATTGTTAGGGCCACCACCTCTGCAGCTGTT AAGGAGGGTGCTCCCTACTTCTGTCCAGGCCTCACCCAGAATCTGACATTCCTGGGATTGTTGTTTCC TCCACATTTGCCCTGCTTGTGTCTGTGTCTAGGTTCTCCCTATTTGTAGGAGCACACATCATCTAAGTG GAGGTCCCTATTCCAAAGGGAAGTCAATTTAGCTAGTGACATCTGCATAGACCTGTCTCCAGATGAAGCC ACACTGAGGCAGGAGACTTCTGGCTTTGGGACACAGTCCCGCTGTAACACGGGGAAGTCAGAGTGTTC GC
OTI Disclaimer:	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in miRBase . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0009965



[View online »](#)