

Product datasheet for **SC400891**

Mir1971 Mouse MicroRNA Expression Plasmid (MI0009970)

Product data:

Product MicroRNA Expression Plasmids

Type:

Species: Mouse

Vector: pCMV-MIR (PCMVMIR)

Precursor: GCAUCAUUUGGUAAGGCUGGGCUGAGAGGACAUCUUUGUCCAUGGACCCAGUAGAUAGGCAGCUUCUGCUUCCAGGAACAACUUGUCUGUUUUUCAGUCUAUC

Mature GUAAAGGCUGGGCUGAGA

Sequence:

Full insert: >MI0009970 GATTCAAATAACTGACAGAGCATACTATTATGTGACATAATAAACCCATTATTGTATAAGGGACTTTC
AGCTGGTATCAAAAGACAGAATTATAAACATTTAGTTATAATCTAATTGGTTTTATTTGTAATTCAT
GGATCTGAGGCAGCCCCATTCTACCAACAGAATGAAAACTCCTGCTACACAATGGCAGAAGAGGGGT
TACAAAGAAATAGGACATGGGAGATAAACCAAGAGGTTGGCATCATGCTACTTCAGGTTAGATGTACATG
CTCTTGATGTACAGGGCATCATTTGGTAAAGGCTGGGCTGAGAGGACATCTTTGTTCCATGGACCCAGTA
GATAGGCAGCTTCTGCTCCAGGAACAACCTTGTCTGTTTTTCAGTCTATCTACTTCCTAAAATTCAATT
TGCTCACATAGCACTTAAGAGAAGAGCTCTCATTTTTGTTTTGCTTTTTGTTGTTCTGCTGTGTTTG
TTTGTCTTCTGTTTTGTTGGGATCTGGTGCAGGAACCTCAGCTGGAAATAATAGGCTTCCATGATTTTTTG
TGTGTTTGTGTTTGTGTTTGTGTTTCAAGACAGGGTCTCACTATGTAGCTCTGACTGTTGGAACCTCAC
TAGGTAGACCAGGTCGGTATTGAACCTAGAGATCTGCTTCTGCTCCAGAGTGCTTGGA

OTI All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence [imiRBase](#).

Disclaimer: The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Kanamycin (25 ug/mL)

Selection:

Mammalian Neomycin

Cell

Selection:

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Both mature and minor (star) miRNA can be over expressed from the construct.

Contains:

RefSeq: [MI0009970](#)



[View online »](#)