

Product datasheet for **SC400899**

Mir19b-1 Mouse MicroRNA Expression Plasmid (MI0000718)

Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Mouse
Vector:	pCMV-MIR (PCMVMIR)
Precursor:	CACUGGUCUAUGGUUAGUUUUGCAGGUUUGCAUCCAGCUGUAUAAUUAUCUGUGUGCAAUCCAUGCAAACUGACUGUGGUGUG
Mature Sequence:	UGUGCAAUCCAUGCAAACUGA
Full insert:	>MI0000718 GCATAGTTGCACTACAAGAAGATGTAGTTGTGCAAATCTATGCAAACCTGATGGTGGCCTGCTATTTAC TTCAAGTGTTGTTTTTTTTAACTAATTTTGTATTTTATTGTGCGATGTAGAGCCTGCGTGGTGTGT GTGATGTGACAGCTTCTGTAGCACTAAAGTGCTTATAGTGCAGGTAGTGTGTAGCCATCTACTGCATTAC GAGCACTTAAAGTACTGCCAGCTGTAGAAGTCCAGCCTCGCCTGGCCATCGCCCAGCCAACTGTCCTGTT ATTGAGCACTGGTCTATGGTTAGTTTTGCAGGTTTGCATCCAGCTGTATAATATTCTGCTGTGCAAATCC ATGCAAACCTGACTGTGGTGGTGAAAAGTCTGTAGAGAAGTAAGGGAAAATCAAACCCCTTCTACACAG GTTGGGATTTGTCGCAATGCTGTGTTTCTCTGTATGGTATTGCACTTGTCCCGCCTGTTGAGTTTGGTG GGGATTGTGACCAGAAGATGTGAAAATGACAATATTGCTGAAGATGCCGATTTCCACTGTAAAAATGTTT ACGAGCTGATGGAGTGTATGAACACATTCTGTAACCTGTGATTCTATTGCTTTCGTTGGTGTGAAAAC GGGATTTTTTAATTGAAATTTGTGTT
OTI Disclaimer:	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in miRBase . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	MI0000568 Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0000718
Synonyms:	Mirn19b-1
Locus ID:	751527



[View online »](#)