

Product datasheet for **SC400911**

Mir208a Mouse MicroRNA Expression Plasmid (MI0000555)

Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Mouse
Vector:	pCMV-MIR (PCMVIR)
Precursor:	UUCCUUUGACGGGUGAGCUUUUUGGCCCGGGUUUACCUGACACUCACGUUAUAGACGAGCAAAAAGCUUGUUGGUCAGAGGAG
Mature Sequence:	AUAAGACGAGCAAAAAGCUUGU
Full insert:	>MI0000555 TACAGCGTGTCTGTCCAAGGCCAACTCAGAGGTGGCCAGTGGAGGACCAAGTATGAGACGGATGCCAT ACAGAGGACGGAGGAGCTGGAGGAAGCCAAGTGAAGTCCACTGTGGCTGCCCTGGATTTGCCAAGAGAAT GGAGCAGAGCGTCTGCTGTGGAGACAGGCAGAGTAGATGCAGGAAAGAGCTTTGGGGTGGGAGAACAGAG GATGGGGAGACCTGGTGGAGCCATGGGAGTGTGGTAGAAAGCCTGAGTCTGCAACCTGTAACGCCTGG CTGGCCCTTGCTTACTTCTTTGACGGGTGAGCTTTTGGCCCGGGTTATACCTGACACTCACGTATAAGA CGAGCAAAAAGCTTGTGGTCAAGGAGCTATCATCATGCAGCCTGGGCTGAGGGAAGCGGGTGGGGT GCACAGTCACACAGATGTGATGTATGTGTGGGAGAGTCTGAGCAGCCTGTGATGGACGTCTCTGTTCTAG GCTCCAGTAGCCCTCACATCCTCTAGGCTGGGGTTCCAACATCTTAGGGAAGTGTGGGAGAAGCTGGGT CTTTCCCATCAACTCGGAAGGAGGAGATCATCTCAGTCTGTGGCTCCAAGTACAGTCACATACAT TCCCCTGGACAGGAAGAAGCTGGCTCAGA
OTI Disclaimer:	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in miRBase . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0000555
Synonyms:	Mirn208 mir-208 Mirn208a mmu-mir-208
Locus ID:	387204



[View online »](#)