

Product datasheet for **SC400916**

Mir210 Mouse MicroRNA Expression Plasmid (MI0000695)

Product data:

Product MicroRNA Expression Plasmids

Type:

Species: Mouse

Vector: pCMV-MIR (PCMVMIR)

Precursor: CCGGGGCAGUCCUCCAGGCUCAGGACAGCCACUGCCCACCGCACACUGCGUUGCUCGACCCACUGUGCGUGACAGCGGCUGAUCUGUCCUGGGCAGCGC

Mature Sequence: CUGUGCGUGUGACAGCGGCUGA

Full insert: >MI0000695 GGGGTAACCTGAACCCGATTGGGGCCCGCCTTGGGTGTGAGTGTGGCCACTGCTGCCACTGAAGTTGATCC
AAGAGCTGGCTATGGCCCTGCGCCGGGGTTCGAGGTGAAATAGAAGGGTTACAAGTTGGAGGATGGGGG
AGCACCTTCGAGGCTGCTTTTCTCAGGGGATATGGGTATTGGGCGGGACAAGAGAGGAGTGGCTCTGTT
GCCAAAGTTAGAGGAGATGAGGGCGGGGGCGGAGAGGAGACCCCTCCACTTGCCACCCCGCGGGG
AGTCCCTCCAGGCTCAGGACAGCCACTGCCACCGCACACTGCGTTGCTCCGGACCCACTGTGCGTGTGA
CAGCGGCTGATCTGTCCCTGGGCAGCGCAACCCCTGGCCCTGTGACCAGAGTCTGGGCGGCGGCAAAT
GGTCCACTCTGGGGCGGAGGAAGGGCACCTGGATTGATGATAGACAGGGTGGTAGCAGGTCAGTGAAGAC
TGTACCTACAGACAGGCCTGCTTGTAGGAGCCAGCTGAGACACGTTGGGCTGGGGGCTTCTCTGTCCTG
CACCTGTGGGTGAGGCTTGGGCCAGATACCTCTTGGAGGGAATTTCTTCTTAGGTGATTAGGGTTC
CAGGCTGTGCAGTCCCAGAAGAGCCACTGTCAAAGTCATCAAGCC

OTI All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence [imiRBase](#). The
Disclaimer: flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms.
This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0000695](#)

Synonyms: Mirn210 mmu-mir-210

Locus ID: 387206



[View online »](#)