

Product datasheet for SC400925

Mir218-2 Mouse MicroRNA Expression Plasmid (MI000701)

Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Mouse
Vector:	pCMV-MIR (PCMVMIR)
Precursor:	GACCAGUUGCCGCGGGCUUCCUUUGUGCUUGAUCUAACCAUGUGGUGGAACGAUGGAAACGGAACAUGGUUCUGUCAAGCACCGCGAAAGCAUCGCUCUCUCU
Mature Sequence:	UUGUGCUUGAUCUAACCAUGU
Full insert:	>MI000701 ACAGGCTCAGTCTGGGTTAGTATGTAATTGGTGTGTCAGCTGCCAGAATGCTCTGGATATGAGTGTAT GAAGGTTTCCCACCTTGTCCCTCTTTGCTGCTGAGTCTCCCTGGGGCAGGGCACTGGGAGTGGGGGCCC CTGTCTTTGGCAGGATCCAGTGGACCGCACTCCTGGTACTAATCACGCTCAGTGGGGCCTGCTCTGGCT TCAGTCTTCCACGCTGCTTCCCTCTGAGCGCTCCTGCTCCTCTCTGACGCTGCTTCCCTGACCTTGACTC TGACCAGTTGCCGCGGGGCTTTCCTTTGTGCTTGATCTAACCATGTGGTGGAACGATGGAACGGAACAT GGTTCTGTCAAGCACCGCGAAAGCATCGCTCTCCTGCAGCATGGCCTGCTGCCACCGCCGCCAC CACCACCACCACCACCGCTGGACTCTTCTCTGCTCTGGATCACCGCGCTGCCAGACCACTTCTCCCAG TCTTAACTCCCAGAGGGAAGGGGGTGGGAATCCTGAAGGAGTCAGGTACATGGTGACAGGAGGTAGAGA AACCCACGGGCGCTGCTTCGAGTATCCACATTACTTAATGAGATGTAGTCACTATGCTCCAAGTCTGGT CTCAAGGGATCCTCAAGAGTAGCAGCTAGGTTCCAATTTC
OTI Disclaimer:	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence imiRBase . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI000701
Synonyms:	Mirn218-2 mir-218-2
Locus ID:	723924



[View online »](#)