

Product datasheet for SC400926

Mir219-1 Mouse MicroRNA Expression Plasmid (MI0000702)

Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Mouse
Vector:	pCMV-MIR (PCMVMIR)
Precursor:	CCGUCCCGGGCCGCGGUCCUGAUUGUCCAAACGCAAUUCUCGAGUCUCUGGCUCGCGCCGAGAGUUGCGUCUGGACGUCCGAGCCGCCCCCAACCCGGA
Mature Sequence:	UGAUUGUCCAAACGCAAUUCU
Full insert:	>MI0000702 CTGGCAGCCGCTGCACAGGCTGGCGGGGGGACGAATACCCATCACGGTCTTCTCCTGTGCCCCCTTAT TCCCGACTCTGAACCGGCTCTCAACCGGTATGGCCCTGCTCCTTCTCCTTCTCTCCTTCTCGTCCGCT TCCTTCACTCCGAGTCCCGGAGTGGCTTTGCCTACTTCCCACGCCAGACATCACTCGTGTGCTCCCCG ACTCCCCCAGGGGCCCGCCCATCTGGTGTCCCGCTTCCGCGGGCGTCTCTATCCCCGGTGCCTCCG ATGCGCTTCTCTCCCGTCCCGGGCCGCGGCTCCTGATTGTCCAACGCAATTCTCGAGTCTCTGGCTCC GGCCGAGAGTTGCGTCTGGACGTCCCGAGCCGCCGCCCAACCTCGAGGGGAGAGGCCGGGCGGAGC GCGGGACGGTCCCCGGACACCCGGGACCCCGCGCAGCCCTAGGATCTCAGGGCGCGCACCCAGCAGCAG AGGTGCACTGCGGGCCGCTTCCGGGGTTCTCAGGCCCGGCTTGGAGAGAAGTTGGGGATCAGCGCCACG GAGGGGATCCTCCGCCCTCGGGCTCTCCACATGCACCAGCACCCACCGGGCCGGCTCACCGCGGCCCGG CGGGGAGGGGGAACCCGGTTCACGGCGCGGTGGCGCGGAGGGCGGTGACGTCTCCGTC
OTI Disclaimer:	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence imiRBase . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0000702
Synonyms:	Mirn219-1
Locus ID:	723823



[View online »](#)