

Product datasheet for **SC400936**

Mir24-1 Mouse MicroRNA Expression Plasmid (MI0000231)

Product data:

Product Type: MicroRNA Expression Plasmids
Species: Mouse
Vector: pCMV-MIR (PCMVIR)
Precursor: CUCCGGUGCCUACUGAGCUGAUUUCAGUUCUCAUUUCACACACUGGCUCAGUUCAGCAGGAACAGGAG
Mature Sequence: GUGCCUACUGAGCUGAUUUCAGU
Full insert: >MI0000231 GACCACCCTTAAGTGTGTCTACAGAAAATACAGTTCTACTGCTAATTAGTTCACTGTAATGTCAAGTGTA
TTGGGAGTCACAGAAAATAAACCAGGTGCATCAAGGAACTGAGCCAACCTTGCTCAGCTGATGCCA
CACGTGATGGGTGGGGTGTGGTTCTTTGTCCCTGCTGGGCGTGGTGAACCTTTCTTGTATTTGCACTC
CAGGTCTCCATGTCTCCACAGTCGCTGGCAAGATGATGCCAGGGCAGCCAGGAGGCACTGCCAGCTG
TGATTGGACCCGCCCTCCGGTGCCTACTGAGCTGATATCAGTTCTCATTTCACACACTGGCTCAGTTTCAG
CAGGAACAGGAGTCGAGCCCTTGAGCAAAAAGCCTACCCATCTGTAAGTGCCTGAGGCAGGCCTCCCACG
GAGAAGAAGATGAAGGGCCTGAGTGTCTGGTACCAGACCAACAGCAGTGGTATCCCTCTGTATGTGCTC
CCGAGGTGGGGATTTGTAGAGTGCAGAGCTAACCTTGGAAAGTTTGCATTCCCTAGAAGCAGGTAGGGCA
GGCGCTGAGCATTGGCCACGTGTTGGGGCAGAGCACACCATGTACCCTGCCTTGAGCAGTGAGAAATGC
CTGTCTGA

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)
Mammalian Cell Selection: Neomycin
Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq: [MI0000231](#)
Synonyms: Mirn24 miR-24 Mirn189 miR-189 Mirn24-1 miR-24-1
Locus ID: 387142



[View online »](#)