

## Product datasheet for **SC400942**

### Mir27a Mouse MicroRNA Expression Plasmid (MI0000578)

#### Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Mouse
Vector:	pCMV-MIR (PCMVMIR)
Precursor:	UGGCCUGAGGAGCAGGGCUUAGCUGCUUGUGAGCAAGGUCCACAGCAAAGUCGUGUUCACAGUGGCUAAGUUCGCCCCUGGACCC
Mature Sequence:	AGGGCUUAGCUGCUUGUGAGCA
Full insert:	>MI0000578 CCTTGTGTTTCAGCTATGTGAGACCCAGCCTGGTCAAGATAGGCAGGCAAGCAAGAATGCTCCAATCTCAC TGTCTCTTCTTCTCTCTTTAGGTGCTACACTCCGCTCCACCTGCTGCCCCCTCACTGTGCCTCGGAC GGCTGGGGTTCCTGGGGATGGGATTTGATGCCAGTCAAAATCACATTGCCAGGGATTTCCAAGTACCC TGTGCTCTGCCTTGGGGGCTCCTGTCGCAAGGATGTCTGTCTTGGGTACTGGGTAACAGAGAAGCCTAT CATGACAACTGGCCTGAGGAGCAGGGCTTAGCTGCTTGTGAGCAAGGTCCACAGCAAAGTCGTGTTTACA GTGGCTAAGTCCGCCCCCTGGACCCATCTCCTCAGGCCGCTGCTGCCACCCCTGCTCTGCCTCCTGGC TGGCAGGGCTCTGCCTCTCTCCGGGCTCCGCTCCCGTGCCTACTGAGCTGAAACAGTTGATTCCAGTGC ACTGGCTCAGTTCAGCAGGAACAGGAGTCCAGCCCCCTAGGAGCTGGCAACCCCTGCAGAATGTCTGCCGT GAGCAAGACCTGGAACAATGGGGTAATTAGTCCAACCTAAGCCCCGTTGTCTGCTGGAGGCTCAGCCTC TAAGCCTGCATCTGGGTGTCCACCTTGAGATTCTGTG
OTI Disclaimer:	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in <a href="#">miRBase</a> . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	<a href="#">MI0000571</a> , <a href="#">MI0000578</a> Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	<a href="#">MI0000578</a>
Synonyms:	Mirn27a mmu-mir-27a
Locus ID:	387220



[View online »](#)