

Product datasheet for **SC400994**

Mir335 Mouse MicroRNA Expression Plasmid (MI0000817)

Product data:

Product Type: MicroRNA Expression Plasmids

Species: Mouse

Vector: pCMV-MIR (PCMVMIR)

Precursor: UCUUUUGGGCGGGGUCAAGAGCAAUAACGAAAAAUGUUUUUUUCGUAAACCGUUUUUCAUUUUGCUCUACCCCUCAUGGGUUUAUGCCA

Mature Sequence: UCAAGAGCAAUAACGAAAAAUGU

Full insert: >MI0000817 TCTTTTGAGCGCCCTAGTGTCAAGTAAGCTATATAGTATTTTCGAGATCTATCTATACCTTTGTTTCCCC
CACGCCCCACCCCATTCCTAAATTCATGCACTATTTCCGTGATTTGTTATATGTAGCTGATAAGATTGA
CTTGATGACTGCTATTTAATTACCACACTCAGGATATTCAGTTATATTTAACCTTCATTATTTTCATTTT
TCACCTAATTTTTACAACAGCATTGGAGATAGGTGCCATTAACCTCATTTTATAGCTGAAGAAACCGAG
AAACAGATTGAAAAGATTTCTTTTGGGCGGGGTCAAGAGCAATAACGAAAAATGTTTGTTCGTAA
ACCGTTTTTCATTATTGCTCCTGACCCCTCTCATGGGTTATAGCCAGTCGAGTAAGTATGCATTTTTTT
TTTTATGTTGGCTTGTGTTCTTTCCAGACTTGATATTCAGTGTTACATCTATTGTCATAGGTTCTTT
AGGTTAAATCCAGTATATGATATTTAATTGCTTTAATTAATAATTTCTGAAAACATGCTACCAAAGT
TTTATCCACTAGCAGTATTAATTTCAATCCAGTACTGTTGGCAAGGATTTTATTGTTTGAATAATCTG
AAATATGTCACATACAGAAAAAGTTCCCATCACTATTAGCCCATCTTG

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0000817](#)

Synonyms: Mirn335 mir-335 mmu-mir-335

Locus ID: 723930



[View online »](#)