

Product datasheet for **SC400998**

Mir340 Mouse MicroRNA Expression Plasmid (MI0000623)

Product data:

Product Type: MicroRNA Expression Plasmids

Species: Mouse

Vector: pCMV-MIR (PCMVMIR)

Precursor: CAUUGUACUUGGUGAUUUAAGCAAUGAGACUGAUUGUCAUUGUCGUUUGUGGGAUCCGUCUCAGUUACUUUJAGCCAUACCGGUAUCUUA

Mature Sequence: UUAUAAAGCAAUGAGACUGAUU

Full insert: >MI0000623 CAGTCTCCTTTGGCAACATTGCCAGTCTCCTCTTTTGGCTAGGATGATAAAGCTACAGTGAGATCTATT
CACTCTTGCTCCTAACAGTGAATATTGACAAGTGTGAACATAATTAATTTTTCTTTCCCAAACAGCT
TCCTGTTGAGATTAGTTGACGCGCCCATTTAGTGTTCGGTTAAGGAAATTACTGGATTGGTTGCCTTGCC
TACCTGTGCAGCTCTCCCACTTCCAGCTTGAGTCTTCAAGAGGGAGCCTCAGAAGTATTCTAAAACC
AATTGTAAGTGGTGTGATTATAAAGCAATGAGACTGATTGTCATATGTCGTTTGTGGGATCCGTCTCAGT
TACTTTATAGCCATACCTGGTATCTTAACACCACAGATCATGCCTGTTGATCAACATTGTAAGGACAAA
ACTTAAGGTTTCTAAATGGAAATTTAAACTGGTTTAGCATTCTTTTATGATTCTTTTAACTCACTCACC
AGAGGCAGATGTTAACTTAATAATGCTAGTTTTAATTTTGTGTCTCCTTTTTGTTTATAGATTTTTTAA
GCAAAAATAAACTGGTAACTACTACTCTAATAAAGTGAATAAACCGTTTTGAATACAGTATTATGCAAA
TATATATTCCATTGAGCTATATTTTTTTCAAACCTCTGTTACT

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0000623](#)

Synonyms: Mirn340 mmu-mir-340

Locus ID: 723845



[View online »](#)