

Product datasheet for **SC401005**

Mir346 Mouse MicroRNA Expression Plasmid (MI0000634)

Product data:

Product Type: MicroRNA Expression Plasmids

Species: Mouse

Vector: pCMV-MIR (PCMVMIR)

Precursor: UCUGUGUUGGGCGUCUGUCUGCCCGAGUGCCUGCCUCUCUGUUGCUCUGAAGGAGGCAGGGCCUGGGCCUGCAGCUGCCUGGGCAGAGCUGCUCCUUC

Mature Sequence: UGUCUGCCCGAGUGCCUGCCUCU

Full insert: >MI0000634 GGATTGTTGAGCAGCCTATGGCTGTCCTAGTGTGTAGCCAAGGCGACATCTTTCTGGTCATGCTACTTT
AGTTACTTAACTAGCCAGTAAAAGCAAAAGGCCAGGCAAGTGCAGGTCCTCCAGGCTTCTGTGAAGGATGGC
ACTACTGGGAGGATTTCTGAGGGTGGCCAAGAAGCAGCTCTTCCCTGCGAGTCTCTGAGTCTCTGCTC
AACAGAGTTGGCCTCCTGGAGAGGCTTCTCATCTGTGATGCCACGGAAGACAGGGCCGAGGTGCTGCGA
CCAGGCGGTTAGGAGGTCTCTGTGTTGGGCGTCTGTCTGCCGAGTGCCTGCCTCTCTGTTGCTCTGAAG
GAGGCAGGGGCTGGGCCTGCAGCTGCCTGGGCAGAGCTGCTCCTTCGAGGGACACACTGCTGGCTCCTCT
CTCTGGTTACAGACTCACGACGTCATGCCTCTGGTCTACAGCACACCTGACCTTACAGCCAGCATGGGC
ATCACTATGCCAAGTCAGGGAAGATCCTGGAAGTCTGGCCCTAAGGACAGAGTGGCGAGGGAGTGGACAT
GACTGAAGATGGGCTCTGAAGAAGTGTAGTGAAGCTTAGGTCCAGTCCACTTCGTCAGATAAAGCTGCAT
CTGGAGTCTTTCATCCAGGGCTGTGTGGCAGGGTAGAGAGCAGCAGGGACTAAA

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0000634](#)

Synonyms: Mirn346 mmu-mir-346

Locus ID: 723847



[View online »](#)