

Product datasheet for SC401009

Mir350 Mouse MicroRNA Expression Plasmid (MI0000640)

Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Mouse
Vector:	pCMV-MIR (PCMVMIR)
Precursor:	AGAUGCCUUGCUCCUACAAGAGUAAAGUGCAUGCGCUUUGGGACAGUGAGGAAAAUAAUGUUCACAAAGCCCAUACACUUUCACCCUUUAGGAGAGUUG
Mature Sequence:	UUCACAAAGCCCAUACACUUUC
Full insert:	>MI0000640 CCCAGCAATAACTCACTTTGTATCCCAAATAGAGCCATGCTATTTATTTCTTCATACTTATTTACAAATG ACATTAACCTTTCTTCAAATTCAGAGCTGTTATCAATAAATAAGTTTCCCATTAGTTTATAAACGGATA CATTAGATTTTGAATAACGGAGGTAGATAGTAGTATAGTAACAATTTATTTAAACTTGGACTTTGAGCAA CTGTGAAAATCAATTGATATGAGATATAGATAGATAGATAGATAGATAGATAGATAGATAGATAGATAGATA TAGATAGATAGATAGATGCCTTCTCTACAAGAGTAAAGTGCATGCGCTTTGGGACAGTGAGGAAAATA ATGTTCAAAAGCCCATACACTTTCCCTTTAGGAGAGTTGCAGCTTGCTTTATACCAACTTCTCACAA AGTACAGCCAAGCTAACCTGAGGGATAGCAAGTGAATGTCTGTATTTAGGATAGTTAGAGATCAGCAGAG GGGTTAAGAAGTATGCTGAGCATTCTGCGTGGTCTTTTACTGAATAGAAACATAAGATTACAACCTCTG ATCTAATCATGCCTCAACTTTCCCTTGGGTAATACACAGAAGCACTAGAAGGACTAATAGCTTTATAA GAAAGAAAATGACTAATCTCCAACAGATTAGTAACAATGAAAC
OTI Disclaimer:	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in miRBase . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0000640
Synonyms:	Mir350 mmu-mir-350
Locus ID:	723921



[View online >](#)