

Product datasheet for **SC401016**

Mir367 Mouse MicroRNA Expression Plasmid (MI0003531)

Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Mouse
Vector:	pCMV-MIR (PCMVIR)
Precursor:	AGGCCGUCACUGUUGCUAACAUCAACUCUGUUGAAUAGAAAUUGGAAUUGCACUUUAGCAAUGGUGAUGGACCU
Mature Sequence:	AAUUGCACUUUAGCAAUGGUGA
Full insert:	>MI0003531 CTCTAAGATGCTGTAGGTTAAGGGAGTTATTGACTGGGCTTCCACCACTTAAACGTGGTTGTAATTGCT TTAGACCTAAGAAAGTAAGTGCTTCCATGTTTTGGTGATGGTATGTTTTCTTTATATTATGCCTTTGT AAAGACCACTCAAGTCTGCTTAAGGGGCTCCTTTACTTTAACATGGAGGCACTTGCTGTGCATTTAAAA ATAAGTGCTTCCATGTTTGAGTGTTGGTGGCTCCTACTTAATCTTCATTGTGCTGCATATAGATAGACAA GGCCGTCAGTGTGCTAACATGCAACTCTGTTGAATAGAAATTGGAATTGCATTTAGCAATGGTGATGG ACCTGTGTCAGTGGCCACTTACTTGGACAGCACTCTCAATGGTATAACCTCCTCGCCATCTTATGTCTC TCCTTGTGTTTCTTCTTATGCTTCTCTTACCTTTGAGAGTGAGTCACCTTTCTGTCTGTCTCCTTCT CCTCCTCAGTAGAGCAGAATTCTAAGCTAGAAAAAGATTTCCGTAAGTCCGTTGTTGTGGTTAAACTG GTTAAACATTTGTAACCCTACTGCCACGGTACCTCTGCTCTCCTTGGAGACCTCGGAAGCTGCTTGGCC GACGCCATCCTTCTGAGCCCTT
OTI Disclaimer:	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in miRBase . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	MI0000402, MI0003531, MI0003718 Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0003531


[View online »](#)