

Product datasheet for **SC401018**

Mir370 Mouse MicroRNA Expression Plasmid (MI0001165)

Product data:

| | |
|---------------------------|---|
| Product Type: | MicroRNA Expression Plasmids |
| Species: | Mouse |
| Vector: | pCMV-MIR (PCMVIR) |
| Precursor: | AGACGGAGAGACCAGGUCACGUCUCUGCAGUUACACAGCUAUGAGUGCCUGCUGGGGUGGAACCGUUUGUCUGUCU |
| Mature Sequence: | GCCUGCUGGGGUGGAACCGGU |
| Full insert: | >MI0001165 GGTAATGTCTCGGCCAAGGTCAGGCATAGAGTTCAGAGCATAGAGCGAGTTTACTGTCTAATTTTTAGG AGACCTGTTGAGCAGATAGGTGTGTGGGTGGGTGGAGCTGTGTGGTATGCAAAGGGTGGGTGTGGCTT TGAGGATCAGTCAGTTTGTCTAGTCTATCCTTCAGTCTGCCTGGTCCATCGTTGGTCCGGTCGGTCCGGTCCG GTAGTCCTTCAGCCTGTTGGACCCAGGAAGGTGCTGCTGTTCTCAGTCGACAAATCCTACGAGTCCGGGG CCTAAGACGGAGAGACCAGGTCACGTCTCTGCAGTTACACAGCTCATGAGTGCCTGCTGGGGTGGAACT GGTTTGTCTGTCTGTCTGCACCATCGCTCGGGTGTGCTGCTGCTGCTGCTGCTGCTGCTGCTGCTACTGCGGGG GGCAAAGATTTGTGAAAGGGGCTTAGAGTTGCTGTGGGAGGGCTTCTGCCTGGCTTTGAGGAGAAGGATT ATCTTTGCTAAGGGATCAGAGGCGTGGATGTTTGGTAGATGGGAACAGAGCACAGAGAGGCTTCTCAATT GATGTAGATGTACCCAGCAAGTTAGAATTAAGGCTGAGTTTGGGGTGGGTATCATTGTGTGGGTATT TTTGTAGGCTATGCTAGCCTCATGACTGTG |
| OTI Disclaimer: | All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in miRBase . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA. |
| E. coli Selection: | Kanamycin (25 ug/mL) |
| Mammalian Cell Selection: | Neomycin |
| Components: | 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube |
| Also Contains: | Both mature and minor (star) miRNA can be over expressed from the construct. |
| RefSeq: | MI0001165 |
| Synonyms: | mmu-mir-370 |
| Locus ID: | 723854 |



[View online »](#)