

## Product datasheet for **SC401021**

### Mir376a Mouse MicroRNA Expression Plasmid (MI0000793)

#### Product data:

**Product Type:** MicroRNA Expression Plasmids  
**Species:** Mouse  
**Vector:** pCMV-MIR (PCMVIR)  
**Precursor:** UAAAAGGUAGAUUCUCCUUCUAUGAGUACAAUUAUUAUGACUAAUCGUAGAGGAAAAUCCACGUUUUC  
**Mature Sequence:** GGUAGAUUCUCCUUCUAUGAGU

**Full insert:** >MI0000793 TCTATGGTTACGTGCTTCTGGATAATCATAGAGGAACATCCACTTTTTTCAGTATCAAATGTCGCTCTGG  
AACTATTTCTGGAAGATTCTGTAGGGTCTGGGGCTGGGGGGAAGATCATGTAATAATGAGTGGGGCGTG  
GCCAGAAATGGCCTTAAGAAGCCAACCTGAGTGTAGTCTGATCTGGGTCTGATATACCGCCTCTTTGGGC  
AGCCTGCTCACTCTGCCAGTGGACTTCCATTTCTGTGACACAAAGCACAAGCTCACAGCATCCTGGATC  
GAGTTGTTCTTTGATATTTAAAAGGTAGATTCTCCTTCTATGAGTACAATATTAATGACTAATCGTAGAG  
GAAAATCCACGTTTTTCAGTATCAAATGCTGCTTGGGAAACCACCGAGGACATGACATCTTACCGAAGAA  
AGCGTGCTATCCGTCGTCTTCTGGGCTTGGCGGGGATAGATAGATTTGTGACTTGTGTGATGTGTGCCT  
ATATGCCAGCCTCTCCAGTGATGTGTTCTGAGACCACTGTGTTCTTTATATCCCTTCTCTGACCTCAA  
TGAAGTTGTTTCTGTGTGCTCACCATGTGTATGTTTTCTGGAGTCGGCTTCTTTTTTCTGTAGCGATA  
TGTCTTGAGAGC

**OTI Disclaimer:** All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

**E. coli Selection:** Kanamycin (25 ug/mL)

**Mammalian Cell Selection:** Neomycin

**Components:** 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

**Also Contains:** Both mature and minor (star) miRNA can be over expressed from the construct.

**RefSeq:** [MI0000793](#)

**Locus ID:** 723855



[View online »](#)